


SI-Aid

Silicone-Faced Wound Dressing



A silicone gel dressing that adheres flexibly and softly and minimizes wound disturbance




- Can be attached quickly and easily and used on various sites by cutting to a suitable shape
- Reduce damages to new tissues with less slippage
- Our silicone adhesion technology minimizes damage to the wound surface and the surrounding skin

Its adhesive properties ensure ease of use.




Since it has enough adherence, you can tape the dressing material without having to hold it in place



- It readily attaches, even to vertical or downward-facing sites, so that anchoring by tape or bandages can be done quickly, even by the patient. It can be used for wounds at virtually any site, since it can be cut to any shape.
- SI-Aid's silicone gel mesh reduces slippage between the dressing and the wound surface, by softly adhering to the skin around the wound. It lessens damage to new tissues and ensures minimum wound disturbance.
- The soft silicone gel mesh disperses the peeling force while pulling away from the skin by contacting the skin in a plane, thus reducing the risk of damage to the wound surface and the surrounding skin.

Instruction for SI-Aid

*Be sure to refer to instruction manual for the use and precautions for use of product. The method of use complies with domestic laws in Japan.



1. Peel off the film

Clean and wash, and disinfect the application site if necessary, then carefully dry the skin surrounding the wound. Hold the tab of the backing film of SI-Aid (white)(blue) and peel it off.



2. Apply

Peel off the backing film without touching the effective surface, and then cover the wound.




3. Anchor

Anchor the SI-Aid with adhesive tape or bandages.



How to peel it off

Slowly peel off the SI-Aid while holding the skin down.



SI-Aid

Type	Product Code No.	Specification (width x length)	Contents Per Box
3	18751	75mm x 75mm	50 pcs
5	18752	100mm x 100mm	50 pcs
7	18753	100mm x 200mm	30 pcs
10	18754	200mm x 300mm	10 pcs
11	18755	300mm x 550mm	5 pcs

EOG Sterilized



ALCARE Co., Ltd.

International Sales Div.

19F, Arcacentral Bldg., 1-2-1 Kinshi, Sumida-ku, Tokyo 130-0013, Japan

TEL.+81-3-5611-7819 FAX.+81-3-5611-7827 MAIL. info@alcare.co.jp

www.alcare.co.jp/e