

# Blood Flow Care

**More safety, easily, and comfortably** —  
**Stagnation of the blood flow in the lower extremities during lying and rest is effectively prevented and the venous reflux is supported.**

Corresponding to various needs and body shapes, the products line provides more appropriate venous return condition.

Medical Stocking  
 (ANSILK-PRO J)



Calf Stocking



Stocking



Calf Stocking(Single Leg)  
 Stocking(Single Leg)

100% Cotton Elastic Bandage  
 (ELASCOT)




## Compression & Gradation Care: Supporting venous return

Compression and Gradation care is a concept of care which ALCARE proposes. It can be achieved by compression and continuous changes of compression gradation.

## Compression Design of ANSILK-PRO J

The three-dimensional structure fitting to the shape of leg and the compression design providing the optimum pressure according to Sigel theory\* (38% increase in blood flow during rest) so it constantly gives adequate compression during lying.




\*Sigel et al, Arch Surg, Vol. 110 Feb. 1975

## Quality of ANSILK-PRO J



Soft touch achieved without using latex  
 Color lines at the ankle and the knee can be a guide for your appropriate application. (Stocking type only)



Two-way delta stretch reduces the feeling of strong pressure on the thigh and assures comfortable fitting. (Stocking type only)



The "comfort heel" increases the smooth feeling at the heel.



A monitor hole enables observation of blood flow conditions.


\*Readjust the position (cover the toe) after observation. Leaving the toe uncovered may result in inadequate blood flow.

## Compression & Gradation Care Products

■ Medical Stocking ■ Elastic Bandage

Stockings and bandages to support venous blood return before, during and after operation.

## Medical Stocking



### ANSILK-PRO J CALF STOCKINGS

Medical Stocking

- The calf stocking type causes no uncomfortable feeling and can be used by many people including those who usually wear no stockings.
- The material with soft texture assures comfortable fitting.
- The graduated compression is based on Sigel Theory.

Type	Product Code No.	Specification		Contents Per Box
		Girth of Ankle	Girth of Calf	
Both Legs	LL 17421	23~27cm	40~46cm	1 pair
	L 17422	21~25cm	36~42cm	
	M 17423	19~23cm	32~38cm	
	S 17424	17~21cm	28~34cm	
Single Leg	SS 17425	15~19cm	24~30cm	1 pce
	LL 17571	23~27cm	40~46cm	
	L 17572	21~25cm	36~42cm	
	M 17573	19~23cm	32~38cm	
	S 17574	17~21cm	28~34cm	
SS 17575	15~19cm	24~30cm		



### ANSILK-PRO J STOCKINGS

Medical Stocking

- It rarely slips down from thigh because of the silicone garter.
- The graduated compression is based on Sigel Theory.
- The two-way delta stretch reduces feeling of strong pressure on the thigh.
- Color lines at the ankle and the knee can be a guide for your appropriate application.

Type	Product Code No.	Specification			Contents Per Box
		Girth of Ankle	Girth of Calf	Girth of Thigh	
Both Legs	LL 17431	23~27cm	40~46cm	55~64cm	1 pair
	L 17432	21~25cm	36~42cm	50~59cm	
	M 17433	19~23cm	32~38cm	47~56cm	
	S 17434	17~21cm	28~34cm	44~53cm	
Single Leg	SS 17435	15~19cm	24~30cm	40~49cm	1 pce
	LL 17581	23~27cm	40~46cm	55~64cm	
	L 17582	21~25cm	36~42cm	50~59cm	
	M 17583	19~23cm	32~38cm	47~56cm	
	S 17584	17~21cm	28~34cm	44~53cm	
SS 17585	15~19cm	24~30cm	40~49cm		

## Elastic Bandage



### ELASCOT

100% Cotton Elastic Bandage

- A 100% cotton elastic bandage has soft touch and excellent absorptivity.
- The original hard twist cotton yarn gives excellent elasticity to the bandage.
- It can be used as a substitute for medical stockings in the case that medical stocking is insufficient.

Type	Product Code No.	Specification		Contents Per Box
		Width	Length (stretched)	
2	11621	5.0cm	4.5m	6 rolls
3	11622	7.5cm	4.5m	6 rolls
4	11623	10.0cm	4.5m	6 rolls
7	11624	17.5cm	4.5m	6 rolls

# ALCARE's "Solution"

for all the people hoping for health

For the medical staff and the patients who seek for Best Medical Care, ALCARE takes a proactive stance on a broad range of issues in various clinical fields.

Pursuing of promotion of patient's early recovery and treatment efficiency

## During / After Operation Management Solution



Hydrocolloid Wound Dressing <KARAYAHESIVE>    Drainage after Operation <SURGIDRAIN ZIPPER>    Pad Attached Moisture Permeable / Waterproof Film Dressing <MULTI FIX-PAD>    Adhesive Dressing for Wound Care <SILKYPORE DRESSING>

Providing optimum materials for prevention and treatment

## Wound Care Solution



Absorbent Pad for Wound <DERMA AID>    Floatation Pad for Decubitus Prevention <BONMAT>

Pursuing a total support of blood flow care


## Compression & Graduation Care Solution



Medical Stocking <ANSILK-II PANTY HOSE>    100% Cotton Elastic Bandage <ELASCOT>    Pad for Compression <CG PAD>    Medical Stocking <ANSILK-PRO J CALF STOCKINGS>

Realizing best care by providing a large variety of products

## Bandage Care Solution



Spandex Elastic Bandage <HIGHSPAN>    Elastic Net Bandage <PRESNET>    Elastic Adhesive Bandage <ELATEX>    Flexible Adhesive Gauze Bandage <SILKYPORE>    Permeable / Waterproof Film Roll <MULTI FIX-ROLL>

Supporting improvement of ostomates' QOL


## Stoma Care Solution



Open-Ended Bag for Colostomy / Ileostomy <LAPACK SUPER G>    Open-Ended Bag for Colostomy / Ileostomy <POSPACK-K>    One-Piece Ostomy System <YOU CARE D>    Two-Piece Ostomy System <PROCARE 2 FA, PROCARE 2 C>

Improving patients' QOL by providing better continence care

## Incontinence Care Solution



Male Incontinence Appliance <YUNIBON>    Male Incontinence Appliance <LEG BAG>

Make it possible to perfectly choose through the stage in prevention, treatment and rehabilitation

## Orthopedic Care Solution



Glass Fiber Casting Tape <CASTLIGHT-α>    Plaster of Paris Bandage <PLASRUN>    Water-Activated Glass Fiber Splint <LIGHT SPLINT-II>    Servical Collar <POLYNECK HARD>



Rib Belt <BUSTBAND-ACE>    Lumbar Supporter <SACRO LIGHT-DX>

Protecting athletes from being injured


## Sports Care Solution



<ANKLE GUARD EX>    <LUMBAR GUARD>    <KNEE GUARD>    <REGUARD SKIN CARE>

## Index

<b>A</b>	
ANSILK-PRO J CALF STOCKINGS.....	6
ANSILK-PRO J STOCKINGS.....	6
<b>C</b>	
COLOCLEAN PC.....	21
COLOCLEAN PC 1-STEP CONE.....	21
<b>D</b>	
DEOFINE POWDER.....	30
DRAINAGE BAG.....	30
<b>E</b>	
ELASCOT.....	6
<b>F</b>	
FLANGE CUTTER.....	28
<b>G</b>	
GASCLEAN.....	30
<b>I</b>	
ILEFINE.....	16
<b>K</b>	
KARAYAHESIVE.....	8
<b>M</b>	
MINI CLOSE.....	26
MINI PAD.....	26
MINI POUCH-K.....	22
MINI POUCH-K SUPER.....	22
MULTI FIX-PAD.....	8
<b>P</b>	
PEDIATRIC PROCARE 1 D.....	17
PEDIATRIC PROCARE 1 POST-OP.....	17
PEDIATRIC PROCARE 1 U.....	24
PEDIATRIC TRUNK BELT.....	29
POSPACK-B.....	13
POSPACK-K.....	13
POSPACK-K STERILE.....	13
POSPACK-LIGHT.....	13
POSPACK-SIMPLE.....	13
PROCARE 1 FREE-CUT C.....	16
PROCARE 1 FREE-CUT D.....	16
PROCARE 1 FREE-CUT U.....	23
PROCARE 1 PRE-CUT C.....	16
PROCARE 1 PRE-CUT D.....	15
PROCARE 1 PRE-CUT Dc.....	16
PROCARE 1 PRE-CUT U.....	23
PROCARE 1 PRE-CUT Uc.....	23
PROCARE 2 BC.....	26
PROCARE 2 C.....	20
PROCARE 2 D.....	20
PROCARE 2 F.....	19
PROCARE 2 FA.....	19
PROCARE 2 FA URO.....	25
PROCARE 2 Fc.....	19
PROCARE 2 Fc URO.....	25
PROCARE 2 IS.....	21
PROCARE 2 M.....	22
PROCARE 2 U.....	25
<b>P</b>	
PROCARE MF PATE.....	27
PROCARE PASTE.....	27
PROCARE POWDER.....	28
PROCARE REMOVER.....	28
PROCARE SOFT WAFER RING.....	27
PROCARE SOFT WAFER STICK.....	27
PROCARE WAFER.....	27
<b>R</b>	
REMOIS CLEANSE.....	28
<b>S</b>	
SD PAD/SD PAD-WP.....	8
SURGIDRAIN OPEN-TOP.....	8
SURGIDRAIN ZIPPER.....	8
<b>T</b>	
TRUNK BELT.....	29
TRUNK BELT WITH FACE PLATE.....	29
<b>U</b>	
URINARY BAG.....	30
<b>Y</b>	
YOU CARE C.....	14
YOU CARE C (10 pieces/box).....	15
YOU CARE Cc.....	15
YOU CARE D.....	14
YOU CARE D (10 pieces/box).....	14
YOU CARE Dc.....	15
YOU CARE TD.....	14
YOU CARE TDc.....	15
YOU CARE U.....	23
YOU CARE Uc.....	23
YOU CARE 2 BC.....	26
YOU CARE 2 Cf.....	19
YOU CARE 2 D.....	19
YOU CARE 2 Df.....	18
YOU CARE 2 F.....	18
YOU CARE 2 Fc.....	18
YOU CARE 2 Fu.....	24
YOU CARE 2 Fuc.....	24
YOU CARE 2 IS.....	21
YOU CARE 2 M.....	22
YOU CARE 2 TDf.....	18
YOU CARE 2 TD.....	18
YOU CARE 2 U.....	24
YOU CARE LEG BAG.....	30
YOU CARE TRUNK BELT.....	29
YOU CARE TRUNK BELT WITH FACE PLATE.....	29



Innovative Solutions for  
Best Care