

# Procare

## Ostomy System

---

Long-selling product  
Since 1991, well-loved by ostomates in the world



## Relief


---

The elastic tape stretches in all direction to accomodate every body movement.

## Comfort

---

The entire film is covered with non-woven material to absorb sweat and obscure excretion.



## Secure

---


"Double-lock structure" is unlikely to cause leakage.

# Procare Ostomy System

## Digestive System Ostomy Appliances


### 1-Piece Ostomy System

#### Procare 1-Pre Cut D ● 3 to 5 Days Exchange




Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
20	14181	20mm φ	14~19mm	300mm×150mm	10 pcs
25	14182	25mm φ	19~24mm	300mm×150mm	10 pcs
30	14183	30mm φ	24~29mm	300mm×150mm	10 pcs
35	14184	35mm φ	29~34mm	300mm×150mm	10 pcs
40	14185	40mm φ	34~39mm	300mm×150mm	10 pcs

#### Procare 1-Pre Cut Dc ● 3 to 5 Days Exchange




Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
22	12366	22mm φ	18~21mm	300mm×150mm	5 pcs
25	12361	25mm φ	21~24mm	300mm×150mm	5 pcs
28	12362	28mm φ	24~27mm	300mm×150mm	5 pcs
32	12363	32mm φ	27~31mm	300mm×150mm	5 pcs
36	12364	36mm φ	31~35mm	300mm×150mm	5 pcs
40	12365	40mm φ	35~39mm	300mm×150mm	5 pcs

#### Procare 1-Free Cut D ● 3 to 5 Days Exchange




Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
50	12331	20mm φ	14~49mm (free-cut)	300mm×150mm	10 pcs

#### Pediatric Procure 1-Post OP ● 1 Day Exchange




Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
—	14211	—	≤34mm (free-cut)	230mm×110mm	10 pcs

#### Procare 1-Pre Cut C ● 1 Day Exchange




Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
20	14191	20mm φ	14~19mm	210mm×150mm	10 pcs
25	14192	25mm φ	19~24mm	210mm×150mm	10 pcs
30	14193	30mm φ	24~29mm	210mm×150mm	10 pcs
35	14194	35mm φ	29~34mm	210mm×150mm	10 pcs
40	14195	40mm φ	34~39mm	210mm×150mm	10 pcs
45	14196	45mm φ	39~44mm	210mm×150mm	10 pcs
50	14197	50mm φ	44~49mm	210mm×150mm	10 pcs

#### Pediatric Procure 1-D ● 2 to 3 Days Exchange



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
—	14221	—	24mm×34mm (oval)	230mm×130mm	10 pcs


#### Procare 1-Free Cut C ● 1 Day Exchange



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
50	12341	20mm φ	14~49mm (free-cut)	210mm×150mm	10 pcs


### 2-Piece Ostomy System

#### Procare 2-F ● 3 to 5 Days Exchange




Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
30	12091	10mm φ	4~29mm (free-cut)	100mm×100mm	30	5 pcs
40	12092	10mm φ	29~39mm (free-cut)	100mm×100mm	40	5 pcs
50	12093	10mm φ	39~49mm (free-cut)	110mm×110mm	50	5 pcs
70	12094	10mm φ	49~69mm (free-cut)	130mm×130mm	70	5 pcs

#### Procare 2-Fc ● 3 to 5 Days Exchange




Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
22-40	12126	22mm φ	18~21mm	110mm φ	40	5 pcs
25-40	12121	25mm φ	21~24mm	110mm φ	40	5 pcs
28-50	12122	28mm φ	24~27mm	120mm φ	50	5 pcs
32-50	12123	32mm φ	27~31mm	120mm φ	50	5 pcs
36-50	12124	36mm φ	31~35mm	120mm φ	50	5 pcs
40-50	12125	40mm φ	35~39mm	120mm φ	50	5 pcs

#### Procare 2-Fa ● 3 to 5 Days Exchange



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
30	12101	10mm φ	4~29mm (free-cut)	100mm φ	30	5 pcs
40	12102	10mm φ	29~39mm (free-cut)	110mm φ	40	5 pcs
50	12103	10mm φ	39~49mm (free-cut)	120mm φ	50	5 pcs

#### Procare 2-D



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
Clear	30	12131	300mm×150mm	10 pcs	
	40	12132	300mm×150mm	10 pcs	
	50	12133	300mm×150mm	10 pcs	
	70	12134	310mm×165mm	10 pcs	
Flesh	30	13641	300mm×150mm	10 pcs	
	40	13642	300mm×150mm	10 pcs	
	50	13643	300mm×150mm	10 pcs	

## Procare 2-C



Type	Product Code No.	Specification		Contents Per Box
		Starter Hole Size	Max Available Size	
Clear	30 12141	10mm $\phi$	4~29mm (free-cut)	10 pcs
	40 12142	10mm $\phi$	29~39mm (free-cut)	10 pcs
	50 12143	10mm $\phi$	39~49mm (free-cut)	10 pcs
Flesh	40 13651	10mm $\phi$	29~39mm (free-cut)	10 pcs
	50 13653	10mm $\phi$	39~49mm (free-cut)	10 pcs

## Procare 2-M



Type	Product Code No.	Specification		Contents Per Box
		Starter Hole Size	Max Available Size	
40	12171	120mm $\phi$	120mm $\phi$	10 pcs
50	12172	120mm $\phi$	120mm $\phi$	10 pcs

## Procare 2-BC



Type	Product Code No.	Specification		Contents Per Box
		Starter Hole Size	Max Available Size	
30	12181	20mm	20mm	5 pcs
40	12182	20mm	20mm	5 pcs
50	12183	20mm	20mm	5 pcs

# Urology System Ostomy Appliances

## 1-Piece Ostomy System

### Procare 1-Pre Cut U 4 to 7 Days Exchange



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
15	14201	15mm $\phi$	9~14mm	260mm×150mm	10 pcs
20	14202	20mm $\phi$	14~19mm	260mm×150mm	10 pcs
25	14203	25mm $\phi$	19~24mm	260mm×150mm	10 pcs
30	14204	30mm $\phi$	24~29mm	260mm×150mm	10 pcs
35	14205	35mm $\phi$	29~34mm	260mm×150mm	10 pcs

### Procare 1-Pre Cut Uc 4 to 7 Days Exchange



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
16	12371	16mm $\phi$	12~15mm	260mm×150mm	5 pcs
19	12372	19mm $\phi$	15~18mm	260mm×150mm	5 pcs
22	12373	22mm $\phi$	18~21mm	260mm×150mm	5 pcs
25	12374	25mm $\phi$	21~24mm	260mm×150mm	5 pcs
28	12375	28mm $\phi$	24~27mm	260mm×150mm	5 pcs

### Procare 1-Free Cut U 4 to 7 Days Exchange



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
50	12351	10mm $\phi$	4~49mm (free-cut)	260mm×150mm	10 pcs

### Pediatric Procure 1-U 2 to 3 Days Exchange



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
—	14231	—	24mm×34mm (oval)	205mm×130mm	10 pcs

## 2-Piece Ostomy System

### Procare 2-Fa Uro 4 to 7 Days Exchange



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
30	13381	10mm $\phi$	4~29mm (free-cut)	100mm $\phi$	30	5 pcs
40	13382	10mm $\phi$	29~39mm (free-cut)	110mm $\phi$	40	5 pcs
50	13383	10mm $\phi$	39~49mm (free-cut)	120mm $\phi$	50	5 pcs

### Procare 2-U



Type	Product Code No.	Specification		Contents Per Box
		Starter Hole Size	Max Available Size	
Clear	30 12151	260mm×150mm	260mm×150mm	10 pcs
	40 12152	260mm×150mm	260mm×150mm	10 pcs
	50 12153	260mm×150mm	260mm×150mm	10 pcs
Flesh	30 13661	260mm×150mm	260mm×150mm	10 pcs
	40 13662	260mm×150mm	260mm×150mm	10 pcs
50 13663	260mm×150mm	260mm×150mm	10 pcs	

### Procare 2-Fc Uro 4 to 7 Days Exchange



Type	Product Code No.	Specification			Contents Per Box	
		Starter Hole Size	Max Available Size	Face Plate Size		
16-40	13391	16mm $\phi$	12~15mm	110mm $\phi$	40	5 pcs
19-40	13392	19mm $\phi$	15~18mm	110mm $\phi$	40	5 pcs
22-40	13393	22mm $\phi$	18~21mm	110mm $\phi$	40	5 pcs
25-40	13394	25mm $\phi$	21~24mm	110mm $\phi$	40	5 pcs
28-50	13395	28mm $\phi$	24~27mm	120mm $\phi$	50	5 pcs



**ALCARE Co., Ltd.**

International Sales Div.  
19F, Arcacentral Bldg., 1-2-1 Kinshi, Sumida-ku, Tokyo  
130-0013, Japan  
TEL.+81-3-5611-7819 FAX.+81-3-5611-7827  
MAIL. [info@alcare.co.jp](mailto:info@alcare.co.jp)  
[www.alcare.co.jp/e](http://www.alcare.co.jp/e)