





















Realization of Active Life for Ostomates


Stoma Care Solution

Since 1965 when "LAPACK" was launched as the first ostomy appliance manufactured in Japan, ALCARE, a stoma care pioneer, has served medical needs and ostomates' demands and made the most of its unique technology to develop and provide various products.

ALCARE continues supplying the optimum solutions to meet a large variety of practical needs in the treatment following ostomy and urinary diversion. ALCARE also provides total stoma care solutions before, during and after operation through social rehabilitation to improve ostomates' QOL and realize their active life.

Information Icon The icons used in this catalogue are listed below. Refer the icon when you find a product with it.

-  **EOG Sterilized** It can be used in an operating room.
-  **Y-ray Sterilized** It can be used in an operating room.
-  **Hypoallergenic Adhesive** The skin friendly hypoallergenic acrylic adhesive
-  **Non-Adherent Pad** The absorptive part of the pad can not adhere to the wound surface.
-  **Water-Proof / Drip-Proof** The water-proof and drip-proof structure prevents contamination from outside.
-  **1-Piece Ostomy System** An ostomy bag and a face plate are integrated into a one-piece appliance.
-  **2-Piece Ostomy System** An ostomy bag and a face plate are separated.
-  **1 Day Exchange**
-  **2 Days Exchange**
-  **2 to 3 Days Exchange** Rough guide for exchange date
-  **3 to 5 Days Exchange**
-  **4 to 7 Days Exchange**
-  **Convex Type** A hard convex plate is incorporated in the face plate to give sufficient height to a stoma.
-  **Safety Plate** A soft plate is incorporated in the face plate to ensure both close contact and flexible fitting to the face plate.
-  **Pre-Cut** A starter hole is cut in advance and scissors are unnecessary when you apply.
-  **Free-Cut** Before applying, cut to fit the stoma size with scissors.
-  **Open-Ended Type** An ostomy bag with a fecal removal opening at the end
-  **Closed-Ended Type** An ostomy bag without a fecal removal opening
-  **Vent for Watery Stool** It is suitable vent for watery to muddy stool.
-  **Vent for Urine** It is suitable vent for urine.



Stoma Care

ALCARE considers improvement of ostomates' QOL as its basic theme of our development and has created many unique products; such as the appliances serving for social rehabilitation and safer skin barrier.

Theme of Development Projects for development of skin barrier that is indispensable for peristomal skin care

ALCARE uses karaya gum as the main raw material for skin barrier to enhance the skin protecting effects.

Features of Karaya Gum

Karaya gum is the natural gum collected from the sap of karaya tree in India. Because of its excellent absorptivity, water-holding capacity, expansivity and safety, karaya gum has been used for various food additives and medicines.


3 Excellent Features

- 1 Buffering Action
- 2 Absorbing Action
- 3 Bacteriostatic Action

Development 1 CMC-Free Skin Barrier

The ingredients of hydrophilic polymer are changed from the combination of karaya gum and CMC to the combination of karaya gum and citrus pectin to improve skin protecting effects.

■ Scanning electron micrographs of the horny layer after one month application of a single hydrophilic polymer (by Tadaoki Kawashima, POLA Chemical Industry Co., Ltd.)




The scanning electron micrographs show the conditions of the horny layers (magnification of 2,000). The conditions of the horny layers treated with karaya gum and citrus pectin were similar to those of the normal horny layer but fissures formed over the horny layers treated with CMC and apple pectin definitely showed aggravation when compared with the conditions of the normal horny layer.

Development 2 Ants' Nest Structure

"Ants' nest structure" of "Procure Barrier for Urinary Care" enables both preservation of skin protection and excellent fitting.


■ Chronological changes of "ants' nest structure" following absorption of water



Development 3 Excellent Buffer Action

"Skin barrier with excellent buffer action" prevents skin irritation resulting from strong alkaline excretion.

The skin barrier with excellent buffer action was developed in collaboration with Toyama Medical and Pharmaceutical University, Japan.



Theme of Development Development of appliances that can be used by the elderly in their self-care processes.


- Products that could be used for a long period of time
- Products that could be used by the elderly in their self-care processes

YOUCARE Series : One-Piece Ostomy System
YOUCARE 2 Series : Two-Piece Ostomy System

Feature 1 Assurance for Skin Friendly

Skin-Friendly
Skin Barrier of Natural Materials

The main ingredients of the skin barrier are the natural materials that rarely irritate the skin such as "karaya gum" and "citrus pectin". Karaya gum keeps the surface of the faceplate weak acid.



Feature 2 Easy and Simple Handling

Secure Fitting without Difficulty

Single Ring Coupling


Aged people having problems in eyes or hands can use YOUCARE 2 Series. The single ring of this series can be fitted smoothly into the groove only by adjusting the two parts and secure fitting is achieved.



User-Friendly

Clip-less Ostomy Appliances

Just turn up the plate along with the blue fold up line and secure the end with Velcro closures attached. Opening and closing are easy. The keep-shaping plate makes discharge and wiping easy.




By folding the keep-shaping plate in V shape, open-end can be kept open.

Feature 3 Comfortable Application

No smell, No leakage

Membrane Charcoal Filter


The new structure consisting of the membrane with excellent air permeability and waterproofing property (air permeable waterproof membrane) and fibrous active carbon with deodorizing effects (charcoal) enables effective protection of accumulated excretion during gas release. Users are no longer afraid of leakage of smell and stool.



Scientific Analysis of Body Movement

Diamond Shape

Anatomically designed Diamond Shape effectively prevents detachment and leakage by increased flexibility following the curved groin and the abdominal movement. Because of the compact design, he/she hardly notices uncomfortable feeling. It offers higher levels of durability and achieves excellent contact.



Stoma Care

Providing a large variety of products based on individual stoma conditions and lifestyles so that ostomates can select optimum ostomy appliances after social rehabilitation.

Selection of Ostomy Appliances after Social Rehabilitation

		Stoma Conditions		
		Stoma With sufficient height Peristomal No wrinkles, no hollows	Stoma Flat or retracted Peristomal Wrinkles and hollows <small>The convex plate inside of the face plate gives sufficient height to the stoma and prevents leakage of excretion.</small>	
Digestive Stoma	One-Piece Ostomy System	Can be changed everyday	 YOU CARE C PROCARE 1 C	 YOU CARE Cc
		Can be used longer days	 YOU CARE TD YOU CARE D PROCARE 1 D ILEFINE D ILEFINE D Cap PEDIATRIC PROCARE 1 D	 YOU CARE TDc YOU CARE Dc PROCARE 1 Dc
	Two-Piece Ostomy System	Can be changed everyday	 YOU CARE 2 F + YOU CARE 2 Cf PROCARE 2 Fa + PROCARE 2 C	 YOU CARE 2 Fc + YOU CARE 2 Cf PROCARE 2 Fc + PROCARE 2 C
		Can be used longer days	 YOU CARE 2 F + YOU CARE 2 Df PROCARE 2 Fa + PROCARE 2 D	 YOU CARE 2 Fc + YOU CARE 2 Df PROCARE 2 Fc + PROCARE 2 D
Urinary Stoma	One-Piece Ostomy System	 YOU CARE U PROCARE 1 U PEDIATRIC PROCARE 1 U	 YOU CARE Uc PROCARE 1 Uc	
	Two-Piece Ostomy System	 YOU CARE 2 Fu + YOU CARE 2 U PROCARE 2 Fa URO + PROCARE 2 U	 YOU CARE 2 Fuc + YOU CARE 2 U PROCARE 2 Fc URO + PROCARE 2 U	

Skin Barrier

For filling the wrinkles and hollows

PROCARE WAFER

PROCARE SOFT WAFER RING

PROCARE SOFT WAFER STICK

Enhancement of adhesion between face plate and skin

PROCARE PASTE

* **Flesh** Flesh type is also available. * Other ostomy appliances are also available.

One-Piece Ostomy System



POSPACK-K STERILE

Open-Ended Bag for Colostomy / Ileostomy

- The EOG sterilized bag is the optimum appliance to be used immediately after operation.
- The open-ended bag with a large opening for removal facilitates easy peristomal treatment.
- Maximum acceptable stoma size up to 109 mm.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
70	15001	20mmφ	14~69mm (free-cut)	400mm×150mm	5 pcs
110	15002	—	≤109mm (free-cut)	400mm×170mm	5 pcs

- EOG Sterilized
- 1-Piece Ostomy System
- 1 Day Exchange
- Free-Cut
- Open-Ended Type



POSPACK-K

Open-Ended Bag for Colostomy / Ileostomy

- The open-ended bag with a large opening for removal facilitates easy peristomal treatment.
- Maximum acceptable stoma size up to 109 mm.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
70	15011	20mmφ	14~69mm (free-cut)	400mm×150mm	10 pcs
110	15012	—	≤109mm (free-cut)	400mm×170mm	10 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Free-Cut
- Open-Ended Type



POSPACK-B

Open-Ended Bag for Colostomy / Ileostomy

- The capacious open-ended bag with a large opening for removal is suitable for those who require greater capacity for accumulation of stool.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
—	15021	20mmφ	14~59mm (free-cut)	400mm×150mm	10 pcs

- 1-Piece Ostomy System
- 3 to 5 Days Exchange
- Free-Cut
- Open-Ended Type



POSPACK-LIGHT

Open-Ended Bag for Colostomy / Ileostomy

- A very thin gentle skin barrier.
- The capacious open-ended bag with a large opening for removal is suitable for those who require greater capacity for accumulation of stool.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
20	15031	20mmφ	14~59mm (free-cut)	400mm×150mm	30 pcs
25	15032	25mmφ	19~24mm	400mm×150mm	30 pcs
30	15033	30mmφ	24~29mm	400mm×150mm	30 pcs
35	15034	35mmφ	29~34mm	400mm×150mm	30 pcs
40	15035	40mmφ	34~39mm	400mm×150mm	30 pcs
45	15036	45mmφ	39~44mm	400mm×150mm	30 pcs
50	15037	50mmφ	44~49mm	400mm×150mm	30 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Pre-Cut
- Free-Cut
- Open-Ended Type



POSPACK-SIMPLE

Open-Ended Bag for Colostomy / Ileostomy

- A very thin gentle skin barrier.
- The capacious open-ended bag with a large opening for removal is suitable for those who require greater capacity for accumulation of stool.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
—	15041	20mmφ	14~49mm (free-cut)	400mm×110mm	30 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Free-Cut
- Open-Ended Type

One-Piece Ostomy System



YOU CARE TD

One-Piece Ostomy System

- A thin and durable skin barrier is used.
- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool. (Without filter type is also available.)
- The clip-less bag can provide easy application and reduce uncomfortable feeling during its application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
Flesh	17441	20mmφ	14~59mm (free-cut)	275mm×150mm	10 pcs
Without Filter	17451	20mmφ	14~59mm (free-cut)	275mm×150mm	10 pcs
20	17461	20mmφ	14~59mm (free-cut)	275mm×150mm	10 pcs
25	17462	25mmφ	19~24mm	275mm×150mm	10 pcs
30	17463	30mmφ	24~29mm	275mm×150mm	10 pcs
35	17464	35mmφ	29~34mm	275mm×150mm	10 pcs
40	17465	40mmφ	34~39mm	275mm×150mm	10 pcs
45	17466	45mmφ	39~44mm	275mm×150mm	10 pcs
50	17467	50mmφ	44~49mm	275mm×150mm	10 pcs

- 1-Piece Ostomy System
- 2 Days Exchange
- Pre-Cut
- Free-Cut
- Open-Ended Type



YOU CARE D

One-Piece Ostomy System

- A thin and durable skin barrier is used.
- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool. (Without filter type is also available.)

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
Flesh	16761	20mmφ	14~59mm (free-cut)	275mm×150mm	30 pcs
Without Filter	16771	20mmφ	14~59mm (free-cut)	275mm×150mm	30 pcs
20	16781	20mmφ	14~59mm (free-cut)	275mm×150mm	30 pcs
25	16782	25mmφ	19~24mm	275mm×150mm	30 pcs
30	16783	30mmφ	24~29mm	275mm×150mm	30 pcs
35	16784	35mmφ	29~34mm	275mm×150mm	30 pcs
40	16785	40mmφ	34~39mm	275mm×150mm	30 pcs
45	16786	45mmφ	39~44mm	275mm×150mm	30 pcs
50	16787	50mmφ	44~49mm	275mm×150mm	30 pcs

- 1-Piece Ostomy System
- 2 Days Exchange
- Pre-Cut
- Free-Cut
- Open-Ended Type



YOU CARE D (10 pieces/box)

One-Piece Ostomy System

- The stoma size can be changed immediately after operation and 10 pieces can be enough at this stage. It is recommended as a starter set.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
Without Filter	16891	20mmφ	14~59mm (free-cut)	275mm×150mm	10 pcs
20	16901	20mmφ	14~59mm (free-cut)	275mm×150mm	10 pcs
25	16902	25mmφ	19~24mm	275mm×150mm	10 pcs
30	16903	30mmφ	24~29mm	275mm×150mm	10 pcs
35	16904	35mmφ	29~34mm	275mm×150mm	10 pcs
40	16905	40mmφ	34~39mm	275mm×150mm	10 pcs
45	16906	45mmφ	39~44mm	275mm×150mm	10 pcs
50	16907	50mmφ	44~49mm	275mm×150mm	10 pcs

- 1-Piece Ostomy System
- 2 Days Exchange
- Pre-Cut
- Free-Cut
- Open-Ended Type



YOU CARE C

One-Piece Ostomy System

- A very thin skin barrier is used.
- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
Flesh	16791	20mmφ	14~59mm (free-cut)	200mm×150mm	30 pcs
20	16801	20mmφ	14~59mm (free-cut)	200mm×150mm	30 pcs
25	16802	25mmφ	19~24mm	200mm×150mm	30 pcs
30	16803	30mmφ	24~29mm	200mm×150mm	30 pcs
35	16804	35mmφ	29~34mm	200mm×150mm	30 pcs
40	16805	40mmφ	34~39mm	200mm×150mm	30 pcs
45	16806	45mmφ	39~44mm	200mm×150mm	30 pcs
50	16807	50mmφ	44~49mm	200mm×150mm	30 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Pre-Cut
- Free-Cut
- Closed-Ended Type

One-Piece Ostomy System



YOU CARE C (10 pieces/box)

One-Piece Ostomy System

- The stoma size can be changed immediately after operation and 10 pieces can be enough at this stage. It is recommended as a starter set.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
20	16911	20mmφ	14~59mm (free-cut)	200mm×150mm	10 pcs
25	16912	25mmφ	19~24mm	200mm×150mm	10 pcs
30	16913	30mmφ	24~29mm	200mm×150mm	10 pcs
35	16914	35mmφ	29~34mm	200mm×150mm	10 pcs
40	16915	40mmφ	34~39mm	200mm×150mm	10 pcs
45	16916	45mmφ	39~44mm	200mm×150mm	10 pcs
50	16917	50mmφ	44~49mm	200mm×150mm	10 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Pre-Cut
- Free-Cut
- Closed-Ended Type



YOU CARE TDC

One-Piece Ostomy System

- A thin and durable skin barrier is used.
- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool.
- The clip-less bag can provide easy application and reduce uncomfortable feeling during its application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
22	17471	22mmφ	18~21mm	275mm×150mm	10 pcs
25	17472	25mmφ	21~24mm	275mm×150mm	10 pcs
28	17473	28mmφ	24~27mm	275mm×150mm	10 pcs
32	17474	32mmφ	27~31mm	275mm×150mm	10 pcs
36	17475	36mmφ	31~35mm	275mm×150mm	10 pcs
40	17476	40mmφ	35~39mm	275mm×150mm	10 pcs
45	17477	45mmφ	39~44mm	275mm×150mm	10 pcs
50	17478	50mmφ	44~49mm	275mm×150mm	10 pcs

- 1-Piece Ostomy System
- 2 Days Exchange
- Convex Type
- Pre-Cut
- Open-Ended Type



YOU CARE Dc

One-Piece Ostomy System

- A thin and durable skin barrier is used.
- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool.
- A convex type face plate is suitable for a flat or a retracted stoma.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
22	16811	22mmφ	18~21mm	275mm×150mm	10 pcs
25	16812	25mmφ	21~24mm	275mm×150mm	10 pcs
28	16813	28mmφ	24~27mm	275mm×150mm	10 pcs
32	16814	32mmφ	27~31mm	275mm×150mm	10 pcs
36	16815	36mmφ	31~35mm	275mm×150mm	10 pcs
40	16816	40mmφ	35~39mm	275mm×150mm	10 pcs
45	16817	45mmφ	39~44mm	275mm×150mm	10 pcs
50	16818	50mmφ	44~49mm	275mm×150mm	10 pcs

- 1-Piece Ostomy System
- 2 Days Exchange
- Convex Type
- Pre-Cut
- Open-Ended Type



YOU CARE Cc

One-Piece Ostomy System

- A very thin skin barrier is used.
- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool.
- A convex type face plate is suitable for a flat or a retracted stoma.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
22	16821	22mmφ	18~21mm	200mm×150mm	10 pcs
25	16822	25mmφ	21~24mm	200mm×150mm	10 pcs
28	16823	28mmφ	24~27mm	200mm×150mm	10 pcs
32	16824	32mmφ	27~31mm	200mm×150mm	10 pcs
36	16825	36mmφ	31~35mm	200mm×150mm	10 pcs
40	16826	40mmφ	35~39mm	200mm×150mm	10 pcs
45	16827	45mmφ	39~44mm	200mm×150mm	10 pcs
50	16828	50mmφ	44~49mm	200mm×150mm	10 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Convex Type
- Pre-Cut
- Closed-Ended Type



PROCARE 1 PRE-CUT D

One-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
20	14181	20mmφ	14~19mm	300mm×150mm	10 pcs
25	14182	25mmφ	19~24mm	300mm×150mm	10 pcs
30	14183	30mmφ	24~29mm	300mm×150mm	10 pcs
35	14184	35mmφ	29~34mm	300mm×150mm	10 pcs
40	14185	40mmφ	34~39mm	300mm×150mm	10 pcs

- 1-Piece Ostomy System
- 3 to 5 Days Exchange
- Pre-Cut
- Open-Ended Type

One-Piece Ostomy System



PROCARE 1 FREE-CUT D

One-Piece Ostomy System

- The soft non-woven fabric is attached to both sides of the ostomy bag.
- Application can be completed by checking the position of the stoma from the slit formed on the front side.
- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
50	12331	20mmφ	14~49mm (free-cut)	300mm×150mm	10 pcs

- 1-Piece Ostomy System
- 3 to 5 Days Exchange
- Free-Cut
- Open-Ended Type



PROCARE 1 PRE-CUT C

One-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
20	14191	20mmφ	14~19mm	210mm×150mm	10 pcs
25	14192	25mmφ	19~24mm	210mm×150mm	10 pcs
30	14193	30mmφ	24~29mm	210mm×150mm	10 pcs
35	14194	35mmφ	29~34mm	210mm×150mm	10 pcs
40	14195	40mmφ	34~39mm	210mm×150mm	10 pcs
45	14196	45mmφ	39~44mm	210mm×150mm	10 pcs
50	14197	50mmφ	44~49mm	210mm×150mm	10 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Pre-Cut
- Closed-Ended Type



PROCARE 1 FREE-CUT C

One-Piece Ostomy System

- The soft non-woven fabric is attached to both sides of the ostomy bag.
- Application can be completed by checking the position of the stoma from the slit formed on the front side.
- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
50	12341	20mmφ	14~49mm (free-cut)	210mm×150mm	10 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Free-Cut
- Closed-Ended Type




PROCARE 1 PRE-CUT Dc

One-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.
- A convex type face plate is suitable for a flat or a retracted stoma.
- TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
22	12366	22mmφ	18~21mm	300mm×150mm	5 pcs
25	12361	25mmφ	21~24mm	300mm×150mm	5 pcs
28	12362	28mmφ	24~27mm	300mm×150mm	5 pcs
32	12363	32mmφ	27~31mm	300mm×150mm	5 pcs
36	12364	36mmφ	31~35mm	300mm×150mm	5 pcs
40	12365	40mmφ	35~39mm	300mm×150mm	5 pcs

- 1-Piece Ostomy System
- 3 to 5 Days Exchange
- Convex Type
- Pre-Cut
- Open-Ended Type



ILEFINE

Ileostomy Bag

- The skin barrier with excellent buffer action prevents skin irritation resulting from strong alkaline excretion.
- The safety plate prevents leakage of watery stool into the face plate.
- *D70: Flat face plate with no safety plate.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
D Cap 30	16991	15mmφ	14~29mm (free-cut)	346mm×148mm	10 pcs
D Cap 40	16992	25mmφ	24~39mm (free-cut)	346mm×148mm	10 pcs
D Cap 50	16993	35mmφ	34~49mm (free-cut)	346mm×148mm	10 pcs
D-30	16981	15mmφ	14~29mm (free-cut)	300mm×148mm	10 pcs
D-40	16982	25mmφ	24~39mm (free-cut)	300mm×148mm	10 pcs
D-50	16983	35mmφ	34~49mm (free-cut)	300mm×148mm	10 pcs
D-70	16984	20mmφ	14~69mm (free-cut)	300mm×148mm	5 pcs

- 1-Piece Ostomy System
- 2 to 3 Days Exchange
- Safety Plate
- Free-Cut
- Open-Ended Type
- Vent for Watery Stool

One-Piece Ostomy System



PEDIATRIC PROCARE 1 POST-OP
One-Piece Ostomy System

- The open-ended bag with a large opening for removal facilitates peristomal treatment.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
—	14211	—	≦34mm (free-cut)	230mm×110mm	10 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Free-Cut
- Open-Ended Type




PEDIATRIC PROCARE 1 D
One-Piece Ostomy System

- The oval face plate can not interfere with the body movement but assures sufficient adhesion.
- PEDIATRIC TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
—	14221	—	24mm×34mm (oval)	230mm×130mm	10 pcs

- 1-Piece Ostomy System
- 2 to 3 Days Exchange
- Free-Cut
- Open-Ended Type

Two-Piece Ostomy System




YOU CARE 2 F
Two-Piece Ostomy System

- A whole skin barrier is recommendable to those who have sensitive skin.
- The thinner margins of the face plate assure close contact with the body and reduce the uncomfortable feeling during application.

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
M22	16591	22mmφ	16~34mm (free-cut)	97mm×105mm	M	5 pcs
M25	16592	25mmφ	21~24mm	97mm×105mm	M	5 pcs
M28	16593	28mmφ	24~27mm	97mm×105mm	M	5 pcs
L32	16594	32mmφ	27~31mm	108mm×117mm	L	5 pcs
L36	16595	36mmφ	31~35mm	108mm×117mm	L	5 pcs
L40	16596	40mmφ	35~39mm	108mm×117mm	L	5 pcs
L45	16597	10mmφ	4~44mm (free-cut)	108mm×117mm	L	5 pcs
65	16598	10mmφ	4~64mm (free-cut)	129mm×140mm	65	5 pcs

- 2-Piece Ostomy System
- 3 to 5 Days Exchange
- Pre-Cut
- Free-Cut



YOU CARE 2 Fc
Two-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.
- A convex type face plate is suitable for a flat or a retracted stoma.

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
M22	16601	22mmφ	18~21mm	110mmφ	M	5 pcs
M25	16602	25mmφ	21~24mm	110mmφ	M	5 pcs
M28	16603	28mmφ	24~27mm	110mmφ	M	5 pcs
L32	16604	32mmφ	27~31mm	120mmφ	L	5 pcs
L36	16605	36mmφ	31~35mm	120mmφ	L	5 pcs
L40	16606	40mmφ	35~39mm	120mmφ	L	5 pcs

- 2-Piece Ostomy System
- 3 to 5 Days Exchange
- Convex Type
- Pre-Cut



YOU CARE 2 Tdf
Two-Piece Ostomy System

- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool.
- The clip-less bag can provide easy application and reduce uncomfortable feeling during its application.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)		
Clear	M	17761	300mm×150mm	10 pcs
	L	17762	300mm×150mm	10 pcs
	65	17763	300mm×150mm	10 pcs
Flesh	M	17771	300mm×150mm	10 pcs
	L	17772	300mm×150mm	10 pcs

- 2-Piece Ostomy System
- Open-Ended Type



YOU CARE 2 TD
Two-Piece Ostomy System

- This system is a drainable bag (open-ended type) without a deodorizing filter.
- The clip-less bag can provide easy application and reduce uncomfortable feeling during its application.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)		
Clear	M	17781	300mm×150mm	10 pcs
	L	17782	300mm×150mm	10 pcs
	65	17783	300mm×150mm	10 pcs
Flesh	M	17791	300mm×150mm	10 pcs
	L	17792	300mm×150mm	10 pcs

- 2-Piece Ostomy System
- Open-Ended Type



YOU CARE 2 Df
Two-Piece Ostomy System

- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)		
Clear	M	16611	300mm×150mm	10 pcs
	L	16612	300mm×150mm	10 pcs
	65	16613	300mm×150mm	10 pcs
Flesh	M	17741	300mm×150mm	10 pcs
	L	17742	300mm×150mm	10 pcs

- 2-Piece Ostomy System
- Open-Ended Type

Two-Piece Ostomy System



YOU CARE 2 D

Two-Piece Ostomy System

- This system is a drainable bag (open-ended type) without a deodorizing filter.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)		
Clear	M	16621	300mm×150mm	10 pcs
	L	16622	300mm×150mm	10 pcs
	65	16623	300mm×150mm	10 pcs
Flesh	M	17751	300mm×150mm	10 pcs
	L	17752	300mm×150mm	10 pcs

- 2-Piece Ostomy System
- Open-Ended Type



YOU CARE 2 Cf

Two-Piece Ostomy System

- The highly functional deodorizing filter enables smooth gas release and prevention of leakage of smell and stool.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)		
Clear	M	16631	210mm×150mm	10 pcs
	L	16632	210mm×150mm	10 pcs
Flesh	M	17801	210mm×150mm	10 pcs
	L	17802	210mm×150mm	10 pcs

- 2-Piece Ostomy System
- Closed-Ended Type



PROCARE 2 F

Two-Piece Ostomy System

- A whole skin barrier type is recommendable to those who have sensitive skin.
- Maximum acceptable stoma size up to 69 mm in diameter.

Type	Product Code No.	Starter Hole Size	Specification			Contents Per Box
			Max Available Size	Face Plate Size	Applicable Bag	
30	12091	10mmφ	4~29mm (free-cut)	100mm×100mm	30	5 pcs
40	12092	10mmφ	29~39mm (free-cut)	100mm×100mm	40	5 pcs
50	12093	10mmφ	39~49mm (free-cut)	110mm×110mm	50	5 pcs
70	12094	10mmφ	49~69mm (free-cut)	130mm×130mm	70	5 pcs

- 2-Piece Ostomy System
- 3 to 5 Days Exchange
- Free-Cut



PROCARE 2 Fa

Two-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.

Type	Product Code No.	Starter Hole Size	Specification			Contents Per Box
			Max Available Size	Face Plate Size	Applicable Bag	
30	12101	10mmφ	4~29mm (free-cut)	100mmφ	30	5 pcs
40	12102	10mmφ	29~39mm (free-cut)	110mmφ	40	5 pcs
50	12103	10mmφ	39~49mm (free-cut)	120mmφ	50	5 pcs

- 2-Piece Ostomy System
- 3 to 5 Days Exchange
- Free-Cut



PROCARE 2 Fc

Two-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.
- A convex type face plate is suitable for a flat or a retracted stoma.
- TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Starter Hole Size	Specification			Contents Per Box
			Max Available Size	Face Plate Size	Applicable Bag	
22-40	12126	22mmφ	18~21mm	110mmφ	40	5 pcs
25-40	12121	25mmφ	21~24mm	110mmφ	40	5 pcs
28-50	12122	28mmφ	24~27mm	120mmφ	50	5 pcs
32-50	12123	32mmφ	27~31mm	120mmφ	50	5 pcs
36-50	12124	36mmφ	31~35mm	120mmφ	50	5 pcs
40-50	12125	40mmφ	35~39mm	120mmφ	50	5 pcs

- 2-Piece Ostomy System
- 3 to 5 Days Exchange
- Convex Type
- Pre-Cut

Two-Piece Ostomy System



PROCARE 2 D

Two-Piece Ostomy System

- Two types are available, a transparent type which allows observation of the stool conditions and a flesh color type which obscures its application.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)		
Clear	30	12131	300mm×150mm	10 pcs
	40	12132	300mm×150mm	10 pcs
	50	12133	300mm×150mm	10 pcs
	70	12134	310mm×165mm	10 pcs
Flesh	30	13641	300mm×150mm	10 pcs
	40	13642	300mm×150mm	10 pcs
	50	13643	300mm×150mm	10 pcs

- 2-Piece Ostomy System
- Open-Ended Type



PROCARE 2 C

Two-Piece Ostomy System

- Two types are available, a transparent type which allows observation of the stool conditions and a flesh color type which obscures its application.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)		
Clear	30	12141	210mm×150mm	10 pcs
	40	12142	210mm×150mm	10 pcs
	50	12143	210mm×150mm	10 pcs
Flesh	30	13651	210mm×150mm	10 pcs
	40	13652	210mm×150mm	10 pcs
	50	13653	210mm×150mm	10 pcs

- 2-Piece Ostomy System
- Closed-Ended Type

ALCARE considers safety and ease of use in the process of designing products and provides a total series of appliances for comfortable colon irrigation.

Irrigation Set



COLOCLEAN PC

Irrigation Set

- This set contains all the appliances needed during and after colon irrigation.
- All parts of the appliances can be put into the attached carrying bag, which is convenient for traveling.

Type	Product Code No.	Specification (Contents of Set)	Contents Per Box
Set	13401	IRRIGATION BAG : 1 pce 2-STEP CONE : 1 pce FACE PLATE : 1 pce DRAINABLE SLEEVE : 5 pcs Carrying bag : 1 pce MINI POUCH-K SUPER : 2 pcs YOU CARE C : 2 pcs	1 set

[Parts]

Product Name	Product Code No.	Specification	Contents Per Box
IRRIGATION BAG	13411	Capacity : 2000ml (with a roll clamp)	1 pce
2-STEP CONE	13412	Silicone (with connector)	1 pce

Product Name	Product Code No.	Specification	Contents Per Box
FACE PLATE	13413	Max Available Size 55mm (with belt)	1 pce
DRAINABLE SLEEVE	13414	Bag Size (L×W) : 900mm×100mm	5 pcs

Stopper for Irrigation




1-STEP CONE

Stopper for Irrigation

- It is recommendable for those who have hard to fit stomas when using a standard 2-STEP CONE.
- A tapering cone suitable for a small stoma.
- It is made of silicon which is flexible and provides a comfortable fit to the stoma.

Type	Product Code No.	Specification	Contents Per Box
—	13415	Silicone (with connector)	1 pce

Sleeve for Irrigation



YOU CARE 2 IS

Two-Piece Ostomy System

- Irrigation can be done without removing any face plate of YOU CARE 2.
- *When irrigation, use the parts of COLOCLEAN PC. :DRAINABLE SLEEVE and 2-STEP CONE or 1-STEP CONE

Type	Product Code No.	Specification Bag Size (L×W)	Contents Per Box
M	16691	900mm×100mm	5 pcs
L	16692	900mm×100mm	5 pcs



PROCARE 2 IS

Two-Piece Ostomy System

- Irrigation can be done without removing any face plate of PROCARE 2.
- *When irrigation, use the parts of COLOCLEAN PC. :DRAINABLE SLEEVE and 2-STEP CONE or 1-STEP CONE

Type	Product Code No.	Specification Bag Size (L×W)	Contents Per Box
40	12161	900mm×100mm	5 pcs
50	12162	900mm×100mm	5 pcs

One-Piece Ostomy System



MINI POUCH-K SUPER

Bag to be used after Irrigation

- It is recommended when the volume of excretion is limited such as after irrigation.
- This small, round-shaped bag offers a snug fit, allowing free body movement.
- Non-woven material applied to the outside of the bag makes the stoma discrete.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size	
—	12881	20mmφ	14~49mm (free-cut)	120mmφ	10 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Free-Cut
- Closed-Ended Type



MINI POUCH-K


Bag to be used after Irrigation

- It is recommended when the volume of excretion is limited such as after irrigation.
- This small, round-shaped bag offers a snug fit, allowing free body movement.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size	
—	12871	20mmφ	14~49mm (free-cut)	120mmφ	30 pcs

- 1-Piece Ostomy System
- 1 Day Exchange
- Free-Cut
- Closed-Ended Type

Two-Piece Ostomy System




YOU CARE 2 M

Two-Piece Ostomy System

- This small, round-shaped bag offers a snug fit, allowing free body movement.
- This system has the highly functional deodorizing filter that enables smooth gas release and prevention of leakage of smell and stool.
- This system can be used in combination with the face plate of YOU CARE 2.

Type	Product Code No.	Specification Bag Size (L×W)		Contents Per Box
		M	L	
M	16671	120mm×130mm		10 pcs
L	16672		120mm×130mm	10 pcs

- 2-Piece Ostomy System
- Closed-Ended Type



PROCARE 2 M

Two-Piece Ostomy System

- This small, round-shaped bag offers a snug fit, allowing free body movement.
- This system has the highly functional deodorizing filter.
- This system can be used in combination with the face plate of PROCARE 2.

Type	Product Code No.	Specification Bag Size		Contents Per Box
		40	50	
40	12171	120mmφ		10 pcs
50	12172		120mmφ	10 pcs

- 2-Piece Ostomy System
- Closed-Ended Type

ALCARE uses the highly durable skin barrier for urinary care and provides the ostomates with a large variety of products so that they can select according to their stoma conditions and lifestyles.

■ One-Piece Ostomy System



YOU CARE U

One-Piece Ostomy System

- The smooth durable thin skin barrier developed specifically for urinary care is used for the face plate.
- The double lock stopper enables simple and easy urinal discharge.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
10	15361	10mmφ	4~59mm (free-cut)	260mm×150mm	10 pcs
20	15362	20mmφ	14~19mm	260mm×150mm	10 pcs
25	15363	25mmφ	19~24mm	260mm×150mm	10 pcs
30	15364	30mmφ	24~29mm	260mm×150mm	10 pcs
35	15365	35mmφ	29~34mm	260mm×150mm	10 pcs
40	15366	40mmφ	34~39mm	260mm×150mm	10 pcs

- 1-Piece Ostomy System
- 2 to 3 Days Exchange
- Pre-Cut
- Free-Cut
- Vent for Urine



YOU CARE Uc

One-Piece Ostomy System

- The smooth durable thin skin barrier developed specifically for urinary care is used for the face plate.
- The double lock stopper enables simple and easy urinal discharge.
- A convex type face plate is suitable for a flat or a retracted stoma.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
13	16401	13mmφ	9~12mm	260mm×150mm	10 pcs
16	16402	16mmφ	12~15mm	260mm×150mm	10 pcs
19	16403	19mmφ	15~18mm	260mm×150mm	10 pcs
22	16404	22mmφ	18~21mm	260mm×150mm	10 pcs
25	16405	25mmφ	21~24mm	260mm×150mm	10 pcs
28	16406	28mmφ	24~27mm	260mm×150mm	10 pcs
32	16407	32mmφ	27~31mm	260mm×150mm	10 pcs

- 1-Piece Ostomy System
- 2 to 3 Days Exchange
- Convex Type
- Pre-Cut
- Vent for Urine



PROCARE 1 PRE-CUT U

One-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.
- The single lock stopper enables simple discharge of urine.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
15	14201	15mmφ	9~14mm	260mm×150mm	10 pcs
20	14202	20mmφ	14~19mm	260mm×150mm	10 pcs
25	14203	25mmφ	19~24mm	260mm×150mm	10 pcs
30	14204	30mmφ	24~29mm	260mm×150mm	10 pcs
35	14205	35mmφ	29~34mm	260mm×150mm	10 pcs

- 1-Piece Ostomy System
- 4 to 7 Days Exchange
- Pre-Cut
- Vent for Urine



PROCARE 1 FREE-CUT U

One-Piece Ostomy System

- The soft non-woven fabric is attached to both sides of the bag.
- Application can be completed by checking the position of the stoma from the slit formed on the front side of the non-woven fabric.
- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.
- The single lock stopper enables simple discharge of urine.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
50	12351	10mmφ	4~49mm (free-cut)	260mm×150mm	10 pcs

- 1-Piece Ostomy System
- 4 to 7 Days Exchange
- Free-Cut
- Vent for Urine



PROCARE 1 PRE-CUT Uc

One-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.
- The single lock stopper enables simple discharge of urine.
- A convex type face plate is suitable for a flat or a retracted stoma.
- TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
16	12371	16mmφ	12~15mm	260mm×150mm	5 pcs
19	12372	19mmφ	15~18mm	260mm×150mm	5 pcs
22	12373	22mmφ	18~21mm	260mm×150mm	5 pcs
25	12374	25mmφ	21~24mm	260mm×150mm	5 pcs
28	12375	28mmφ	24~27mm	260mm×150mm	5 pcs

- 1-Piece Ostomy System
- 4 to 7 Days Exchange
- Convex Type
- Pre-Cut
- Vent for Urine

■ One-Piece Ostomy System



PEDIATRIC PROCARE 1 U


One-Piece Ostomy System

- The oval face plate can not interfere with the body movement but assures sufficient adhesion.
- PEDIATRIC TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
—	14231	—	24mm×34mm (oval)	205mm×130mm	10 pcs

- 1-Piece Ostomy System
- 2 to 3 Days Exchange
- Free-Cut
- Vent for Urine

■ Two-Piece Ostomy System




YOU CARE 2 Fu

Two-Piece Ostomy System

- A whole skin barrier type is recommendable to those who have sensitive skin.
- The thinner margins of the face plate assure close contact with the body and reduce the uncomfortable feeling during application.

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
M13	16641	13mmφ	7~34mm (free-cut)	97mm×105mm	M	5 pcs
M16	16642	16mmφ	12~15mm	97mm×105mm	M	5 pcs
M19	16643	19mmφ	15~18mm	97mm×105mm	M	5 pcs
M22	16644	22mmφ	18~21mm	97mm×105mm	M	5 pcs
M25	16645	25mmφ	21~24mm	97mm×105mm	M	5 pcs
M28	16646	28mmφ	24~27mm	97mm×105mm	M	5 pcs
L32	16647	32mmφ	27~44mm (free-cut)	108mm×117mm	L	5 pcs

- 2-Piece Ostomy System
- 4 to 7 Days Exchange
- Pre-Cut
- Free-Cut




YOU CARE 2 Fuc

Two-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.
- A convex type face plate is suitable for a flat or a retracted stoma.

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
M13	16651	13mmφ	9~12mm	110mmφ	M	5 pcs
M16	16652	16mmφ	12~15mm	110mmφ	M	5 pcs
M19	16653	19mmφ	15~18mm	110mmφ	M	5 pcs
M22	16654	22mmφ	18~21mm	110mmφ	M	5 pcs
M25	16655	25mmφ	21~24mm	110mmφ	M	5 pcs
M28	16656	28mmφ	24~27mm	110mmφ	M	5 pcs
L32	16657	32mmφ	27~31mm	120mmφ	L	5 pcs

- 2-Piece Ostomy System
- 4 to 7 Days Exchange
- Convex Type
- Pre-Cut



YOU CARE 2 U

Two-Piece Ostomy System

- The double lock stopper enables simple and easy urinal discharge.
- YOU CARE TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification		Contents Per Box
		Starter Hole Size	Bag Size (L×W)	
M	16661	—	260mm×150mm	10 pcs
L	16662	—	260mm×150mm	10 pcs

- 2-Piece Ostomy System
- Vent for Urine

Urology System Ostomy Appliances

Two-Piece Ostomy System

Two-Piece Ostomy System



PROCARE 2 Fa URO

Two-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	
30	13381	10mmφ	4~29mm (free-cut)	100mmφ	30
40	13382	10mmφ	29~39mm (free-cut)	110mmφ	40
50	13383	10mmφ	39~49mm (free-cut)	120mmφ	50

- 2-Piece Ostomy System
- 4 to 7 Days Exchange
- Free-Cut



PROCARE 2 Fc URO

Two-Piece Ostomy System

- The tape around the face plate prevents maceration and collapse of the margins of the skin barrier during bathing.
- A convex type face plate is suitable for a flat or a retracted stoma.
- TRUNK BELT, P29, is available for more secure application.

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	
16-40	13391	16mmφ	12~15mm	110mmφ	40
19-40	13392	19mmφ	15~18mm	110mmφ	40
22-40	13393	22mmφ	18~21mm	110mmφ	40
25-40	13394	25mmφ	21~24mm	110mmφ	40
28-50	13395	28mmφ	24~27mm	120mmφ	50

- 2-Piece Ostomy System
- 4 to 7 Days Exchange
- Convex Type
- Pre-Cut



PROCARE 2 U

Two-Piece Ostomy System

- Two types are available, a transparent type which allows observation of the urine conditions and a flesh color type which obscures its application.
- The single lock stopper enables simple discharge of urine.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)	Starter Hole Size	
Clear	30	260mm×150mm	10mmφ	10 pcs
	40	260mm×150mm	10mmφ	10 pcs
	50	260mm×150mm	10mmφ	10 pcs
Flesh	30	260mm×150mm	10mmφ	10 pcs
	40	260mm×150mm	10mmφ	10 pcs
	50	260mm×150mm	10mmφ	10 pcs


- 2-Piece Ostomy System
- Vent for Urine

Appliances for Bathing

One-Piece Ostomy System ■ Two-Piece Ostomy System ■ Plaster with Pad

The compact and indistinctive appliances effectively support the ostomates while they enjoy bathing and swimming.

One-Piece Ostomy System



MINI CLOSE

Bag for Bathing

- Use of this type of bag is recommended when the volume of excretion is limited such as after irrigation and during bathing.
- A hypo-allergenic adhesive is used for the face plate.
- This small, oval-shaped bag offers a snug fit, allowing free body movement.
- The flesh color film makes the stoma discreet.

- 1-Piece Ostomy System
- Pre-Cut
- Closed-Ended Type

Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L×W)	
M	12641	40mmφ	≤39mm	120mm×110mm	30 pcs
L	12642	40mmφ	≤39mm	140mm×130mm	30 pcs

Two-Piece Ostomy System



YOU CARE 2 BC

Two-Piece Ostomy System

- A compact dome shape suitable for bathing and swimming.
- This plastic dome is highly durable and allows you to wipe off water on the surface.
- The flesh color film makes the stoma discreet.
- This system can be used in combination with the face plate of YOU CARE 2.

- 2-Piece Ostomy System
- Closed-Ended Type

Type	Product Code No.	Specification (Height of Dome)	Contents Per Box
M	16681	20mm	2 pcs
L	16682	20mm	2 pcs
65	16683	20mm	2 pcs



PROCARE 2 BC

Two-Piece Ostomy System

- A compact dome shape suitable for bathing and swimming.
- This plastic dome is highly durable and allows you to wipe off water on the surface.
- The flesh color film makes the stoma discreet.
- This system can be used in combination with the face plate of PROCARE 2.

- 2-Piece Ostomy System
- Closed-Ended Type

Type	Product Code No.	Specification (Height of Dome)	Contents Per Box
30	12181	20mm	5 pcs
40	12182	20mm	5 pcs
50	12183	20mm	5 pcs

Plaster with Pad



MINI PAD

Pad for Bathing

- Use of this pad is recommended during bathing or at times when transient protection of the stoma is needed.
- An absorbing pad offering friendly protection of the stoma.
- The flesh color base material makes the stoma discreet.

Type	Product Code No.	Specification		Contents Per Box
		Size of Adhesive Part	Size of Absorbing Part	
—	12651	100mm×100mm	65mm×50mm	30 pcs

Accessories

■Skin Barrier

ALCARE provides the skin barrier containing the natural raw materials with excellent skin protecting and a large variety of pouch accessories so that the ostomates can enjoy more comfortable daily life.

■Skin Barrier




PROCARE WAFER

Wafer Skin Barrier (Karaya Gum + Citrus Pectin)

- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- It can be freely cut to the shape of the hollow or the wrinkled area, or it can be used in layers if necessary.

Type	Product Code No.	Specification (L×W)	Contents Per Box
10	12581	100mm×100mm	10 pcs
20	12582	200mm×195mm	5 pcs



PROCARE SOFT WAFER RING

Wafer Skin Barrier (Karaya Gum + Citrus Pectin)

- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- A moderately soft wafer that can be molded with the hands (no scissors needed).
- A ring-shaped wafer that is easy to apply to the peristomal skin area.
- A double-sides adhesive wafer that adheres well to the face plate of various appliances.

Type	Product Code No.	Specification (Diameter of Hole)	Contents Per Box	
Ring	20	15591	20mmφ	8 pcs
	30	15592	30mmφ	8 pcs



PROCARE SOFT WAFER STICK

Wafer Skin Barrier (Karaya Gum + Citrus Pectin)

- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- An extremely soft wafer that can be freely molded.
- A double-sides adhesive wafer that adheres well to the face plate of various appliances.
- Ideal for application to uneven hard-to-apply surfaces such as the area around the fistula.

Type	Product Code No.	Specification (L×W)	Contents Per Box
Stick	15601	15mm×100mm	10 pcs



PROCARE MF PATE

Paste Skin Barrier (Karaya Gum + Citrus Pectin)

- It fills the gap between the stoma and face plate, protecting the skin from discharge.
- An alcohol-free paste which will not make the skin sting but will remain soft.

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	12601	100g	1 tube



PROCARE PASTE

Paste Skin Barrier (Karaya Gum + Citrus Pectin)

- By applying this paste on the inside of the face plate, the adhesion of appliances can be improved.
- A laminate tube keeps the paste moist, making it easy to squeeze the tube down to the last.

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	14431	50g	1 tube

Accessories

■Skin Barrier ■Stoma Care Accessory

■Skin Barrier



PROCARE POWDER

Powder Skin Barrier (Karaya Gum + Citrus Pectin)

- This powder fills the gap between the stoma and the face plate, protecting the skin from discharge.
- It is easy to apply, even to a small gap. The powder offers excellent water absorption.

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	12611	50g	1 bottle

■Stoma Care Accessory



REMOIS CLEANSE

Skin Moisture Cleansing Cream

- Natural oil-in makes the dirt detach from skin and just wiping off makes the skin clean.
- Cleansing with this cream assures refreshed feeling and ostomy appliances can be applied immediately after wiping off the cream properly.
- Moisturizer added cream will make the skin moist and prevent from dryness.

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
Regular	17701	180g	1 pce
Handy	17702	5g / pack	10 pks



PROCARE REMOVER

Ostomy Appliance Remover

- It reduces the peeling pain when removing the face plate of ostomy appliance from the skin.
- It can be effectively used when an appliance needs immediate removal.
- It is individually packed in a foil pack.

Type	Product Code No.	Specification (Fluid Measure)	Contents Per Box
—	13871	2cc	50 pcs



FLANGE CUTTER

Cutter for Flange

- A starter hole on the flange can be easily made using this FLANGE CUTTER.
- The finishes, smooth cut edges of the flange reduce aggravation of the stoma.

Type	Product Code No.	Specification (Length)	Contents Per Box
—	12661	Main Body 123mm	20 pcs

Accessories

Stoma Care Accessory

Stoma Care Accessory



TRUNK BELT WITH FACE PLATE

Belt for One-Piece Bag

- Provides a securer application, preventing the flange from slipping or falling off with the body movement.
- It is recommendable for active athletes or for those who have difficulty in adhering the appliances to the skin.

Type	Product Code No.	Specification		Contents Per Pack
		Waist Size	Applicable Bag	
—	12671	60~110cm (diameter of hole 63mmφ)	PROCARE 1 D,C,U	1 set



YOUCARE TRUNK BELT WITH FACE PLATE

Belt for One-Piece Bag

- Provides a securer application, preventing the flange from slipping or falling off with the body movement.
- It is recommendable for active athletes or for those who have difficulty in adhering the appliances to the skin.

Type	Product Code No.	Specification		Contents Per Pack
		Waist Size	Applicable Bag	
—	15081	60~110cm (diameter of hole 73mmφ)	YOU CARE TD YOU CARE D YOU CARE C YOU CARE U POSPACK-B POSPACK-LIGHT	1 set



TRUNK BELT

Belt for Ostomy Appliance

- Provides a securer application, preventing the flange from slipping or falling off with the body movement.
- Most suitable for active athletes or for those who have difficulty in adhering the appliances to the skin.

Type	Product Code No.	Specification		Contents Per Pack
		Waist Size	Applicable Bag	
1	12691	60~110cm	PROCARE 2 Fc PROCARE 2 Fc URO PROCARE 1 Dc PROCARE 1 Uc	1 pce



YOUCARE TRUNK BELT

Belt for Ostomy Appliance

- Provides a securer application, preventing the flange from slipping or falling off with the body movement.
- Most suitable for active athletes or for those who have difficulty in adhering the appliances to the skin.
- *Not applicable to ILEFINE D-70

Type	Product Code No.	Specification		Contents Per Pack
		Waist Size	Applicable Bag	
—	15091	60~110cm	YOU CARE 2 TDf YOU CARE 2 TD YOU CARE 2 Df YOU CARE 2 D YOU CARE 2 Cf YOU CARE 2 U YOU CARE 2 IS YOU CARE 2 M YOU CARE 2 BC YOU CARE TDc YOU CARE Dc YOU CARE Cc YOU CARE Uc ILEFINE D CAP ILEFINE D	1 pce



PEDIATRIC TRUNK BELT

Belt for One-Piece Bag

- This belt can be used with PEDIATRIC PROCARE 1 only.
- Provides a securer application, preventing the flange from slipping or falling off with the body movement.
- The hooks at both ends provide easy handling. They are to be set at the belt holders for stabilization and are to be removed for detachment.

Type	Product Code No.	Specification		Contents Per Pack
		Waist Size	Applicable Bag	
—	14441	30~50cm	PEDIATRIC PROCARE 1 D PEDIATRIC PROCARE 1 U	1 set

Accessories

Stoma Care Accessory

Stoma Care Accessory




DEOFINE POWDER

Powdery Deodorant For Ostomates

- Plant ingredient inhibits the generation and diffusion of odors.
- Refreshing mint ingredient makes stoma care more comfortable.
- Stick-type individual packaging can provide portability.

Type	Product Code No.	Specification (Net Weight)		Contents Per Box
		1.5g / pack		
—	17891			30 pcs



GASCLEAN

Deodorant Gas Filter

- Prevents the bag from swelling with gas and the pressure on the adhesive part is effectively reduced.
- A charcoal filter deodorizes the released gas.
- Self-adhesive makes for easy handling.

Type	Product Code No.	Specification		Contents Per Box
		Outer Diameter	Thickness	
—	12681	30mmφ	3mm	120 pcs



DRAINAGE BAG

Bag for Drainage

- This bag can be connected to SURGIDRAIN Series.
- The bag can be used during sleep by connecting it to ILEFINE D Cap of the bag for watery to soft stool.
- *Concomitant use of the drainage bag is recommended in case of watery stool. Its concomitant use in case of muddy to soft stool may block the tube and prevent smooth flow.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)	Capacity	
—	15841	320mm×230mm (tube length : 1m)	2000ml	5 pcs



YOUCARE LEG BAG

Leg Bag for Urinary Bag

- A leg bag that connects with the ostomy bag for an easy and convenient way to increase the capacity of urine collection.
- A leg bag that collects urine at the lower leg and reduces weight on the bag adhesive.
- A non-woven lining on the skin-side of the leg bag absorbs sweat and offers greater wearer comfort.

Type	Product Code No.	Specification (Capacity)		Contents Per Box
		500ml		
—	15371			5 pcs



URINARY BAG

Closed-Ended Urinary Bag

- A urinary bag that connects to all types of stoma bags for easy and convenient collection of urine during the night.
- A one-way valve prevents back-flow of urine into the bag.
- An easy, one-touch connection mechanism with ALCARE's urinary pouch.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (L×W)	Capacity	
—	12071	320mm×230mm (tube length : 1m)	2000ml	5 pcs

ALCARE's "Solution"

for all the people hoping for health

For the medical staff and the patients who seek for Best Medical Care, ALCARE takes a proactive stance on a broad range of issues in various clinical fields.

Pursuing of promotion of patient's early recovery and treatment efficiency

During / After Operation Management Solution



Hydrocolloid Wound Dressing <KARAYAHESIVE> Drainage after Operation <SURGIDRAIN ZIPPER> Pad Attached Moisture Permeable / Waterproof Film Dressing <MULTI FIX-PAD> Adhesive Dressing for Wound Care <SILKYPORE DRESSING>

Providing optimum materials for prevention and treatment

Wound Care Solution



Absorbent Pad for Wound <DERMA AID> Floatation Pad for Decubitus Prevention <BONMAT>

Pursuing a total support of blood flow care

Compression & Graduation Care Solution



Medical Stocking <ANSILK-II PANTY HOSE> 100% Cotton Elastic Bandage <ELASCOT> Pad for Compression <CG PAD> Medical Stocking <ANSILK-PRO J CALF STOCKINGS>

Realizing best care by providing a large variety of products

Bandage Care Solution



Spandex Elastic Bandage <HIGHSPAN> Elastic Net Bandage <PRESNET> Elastic Adhesive Bandage <ELATEX> Flexible Adhesive Gauze Bandage <SILKYPORE> Permeable / Waterproof Film Roll <MULTI FIX-ROLL>

Supporting improvement of ostomates' QOL


Stoma Care Solution



Open-Ended Bag for Colostomy / Ileostomy <LAPACK SUPER G> Open-Ended Bag for Colostomy / Ileostomy <POSPACK-K> One-Piece Ostomy System <YOU CARE D> Two-Piece Ostomy System <PROCARE 2 FA, PROCARE 2 C>

Improving patients' QOL by providing better continence care

Incontinence Care Solution



Male Incontinence Appliance <YUNIBON> Male Incontinence Appliance <LEG BAG>

Make it possible to perfectly choose through the stage in prevention, treatment and rehabilitation

Orthopedic Care Solution



Glass Fiber Casting Tape <CASTLIGHT-α> Plaster of Paris Bandage <PLASRUN> Water-Activated Glass Fiber Splint <LIGHT SPLINT-II> Servical Collar <POLYNECK HARD>



Rib Belt <BUSTBAND-ACE> Lumbar Supporter <SACRO LIGHT-DX>

Protecting athletes from being injured

Sports Care Solution



<ANKLE GUARD EX> <LUMBAR GUARD> <KNEE GUARD> <REGUARD SKIN CARE>

Index

A	
ANSILK-PRO J CALF STOCKINGS	6
ANSILK-PRO J STOCKINGS	6
C	
COLOCLEAN PC	21
COLOCLEAN PC 1-STEP CONE	21
D	
DEOFINE POWDER	30
DRAINAGE BAG	30
E	
ELASCOT	6
F	
FLANGE CUTTER	28
G	
GASCLEAN	30
I	
ILEFINE	16
K	
KARAYAHESIVE	8
M	
MINI CLOSE	26
MINI PAD	26
MINI POUCH-K	22
MINI POUCH-K SUPER	22
MULTI FIX-PAD	8
P	
PEDIATRIC PROCARE 1 D	17
PEDIATRIC PROCARE 1 POST-OP	17
PEDIATRIC PROCARE 1 U	24
PEDIATRIC TRUNK BELT	29
POSPACK-B	13
POSPACK-K	13
POSPACK-K STERILE	13
POSPACK-LIGHT	13
POSPACK-SIMPLE	13
PROCARE 1 FREE-CUT C	16
PROCARE 1 FREE-CUT D	16
PROCARE 1 FREE-CUT U	23
PROCARE 1 PRE-CUT C	16
PROCARE 1 PRE-CUT D	15
PROCARE 1 PRE-CUT Dc	16
PROCARE 1 PRE-CUT U	23
PROCARE 1 PRE-CUT Uc	23
PROCARE 2 BC	26
PROCARE 2 C	20
PROCARE 2 D	20
PROCARE 2 F	19
PROCARE 2 FA	19
PROCARE 2 FA URO	25
PROCARE 2 Fc	19
PROCARE 2 Fc URO	25
PROCARE 2 IS	21
PROCARE 2 M	22
PROCARE 2 U	25
P	
PROCARE MF PATE	27
PROCARE PASTE	27
PROCARE POWDER	28
PROCARE REMOVER	28
PROCARE SOFT WAFER RING	27
PROCARE SOFT WAFER STICK	27
PROCARE WAFER	27
R	
REMOIS CLEANSE	28
S	
SD PAD/SD PAD-WP	8
SURGIDRAIN OPEN-TOP	8
SURGIDRAIN ZIPPER	8
T	
TRUNK BELT	29
TRUNK BELT WITH FACE PLATE	29
U	
URINARY BAG	30
Y	
YOU CARE C	14
YOU CARE C (10 pieces/box)	15
YOU CARE Cc	15
YOU CARE D	14
YOU CARE D (10 pieces/box)	14
YOU CARE Dc	15
YOU CARE TD	14
YOU CARE TDc	15
YOU CARE U	23
YOU CARE Uc	23
YOU CARE 2 BC	26
YOU CARE 2 Cf	19
YOU CARE 2 D	19
YOU CARE 2 Df	18
YOU CARE 2 F	18
YOU CARE 2 Fc	18
YOU CARE 2 Fu	24
YOU CARE 2 Fuc	24
YOU CARE 2 IS	21
YOU CARE 2 M	22
YOU CARE 2 TDf	18
YOU CARE 2 TD	18
YOU CARE 2 U	24
YOU CARE LEG BAG	30
YOU CARE TRUNK BELT	29
YOU CARE TRUNK BELT WITH FACE PLATE	29



Innovative Solutions for
Best Care