

Quickfix-N

Catheter Fixation Assisted Tape

Features

- As shown on the guides on the release paper, it is an integrated tape that can be fixated firmly with one "Fixation" and the catheter position can be adjusted at the "Adjustment" position.
- Adhesive tape with low adhesive residue is used to fixate the catheter on one side, making it easy to remove and facilitate daily observation and replacement.



Type	Product Code No.	Specification		Contents Per Box
		Size Width x Length	Applicable Catheter Size	
N2	20361	2.5 cm × 8.5 cm	18 Fr or less	100 pieces