

Pursuing of promotion of patients'
early recovery and treatment efficiency
by supporting from the body surface

During/After Operation Management Solution

During / After Operation Management Solution

Blood Flow Care
Wound Care
Drainage Care

ALCARE

ALCARE Co., Ltd.
INTERNATIONAL SALES DIVISION

19F, Arcacentral Bldg.
1-2-1 Kinshi, Sumida-ku
Tokyo 130-0013, Japan
TEL:81-3-5611-7819
FAX:81-3-5611-7827
<http://www.alcare.co.jp/e>

Printed in Japan
Aug. 2005

ALCARE

Innovative Solutions for
Best Care

Providing best care for all the people hoping for health Innovative Solutions for Best Care

ALCARE bridges between medical care and people, actively addresses the challenges existing between them and creates new values.

The half-century history of ALCARE is a series of technological progress and challenge of product development based on the fundamental principle of customer loyalty.

ALCARE now attempts to reach the next stage with the concept of truly affluent medical and welfare society leading to a better future in mind.

Under the corporate message, "Innovative Solutions for Best Care", ALCARE proposes really innovative solutions to create best care for all the people that provide medical care, need medical care and hope for promotion of health.



ALCARE is actively involved in solutions in the following fields:

- Pursuing of promotion of patients' early recovery and treatment efficiency
- **During / After Operation Management Solution**
- Wound Care Solution
- Compression & Gradation Care Solution
- Bandage Care Solution
- Stoma Care Solution
- Orthopedic Care Solution
- Sports Care Solution

ALCARE keeps challenging to achieve our dream.

During / After Operation Management Solution


C O N T E N T S

Introduction of During / After Operation Management Solution.....	3
■ [Blood Flow Care]	5
● Compression & Gradation Care Products	
Medical Stocking.....	6
Elastic Bandage.....	6
■ [Wound Care]	7
● Wound Care Products	
Hydrocolloid Dressing.....	9
Adhesive Island Dressing.....	9
Absorbent Pad.....	10
Film Dressing.....	10
Fixation Material.....	10
■ [Drainage Care]	11
● Drainage Care Products	
Drainage Pouch.....	13
Drainage Pad.....	13
Absorbent Pad.....	13
Drainage Bag.....	14
Accessory.....	14
Skin Barrier for Post-PEG Treatment.....	14
Products Introduction by Solution.....	15
Index.....	17

Safer and simpler care during operation for patients' optimum comfort

During/After Operation Management Solution

ALCARE has introduced new product concepts to develop and provide the alternatives to gauze as an important medical material that plays a central role in wound and drainage care. Wound care is changing from treatment with gauze to wet dressing. Drainage care is changing from treatment with gauze to pouch drainage. With promotion of patients' early recovery and treatment efficiency in during / after-operative care in mind, ALCARE provides value-added solutions in conjunction with blood flow care and urination care.



- Information Icon** The icons used in this catalogue are listed below. Refer the icon when you find a product with it.
- EOG Sterilized** It can be used in an operating room.
 - γ-ray Sterilized** It can be used in an operating room.
 - Hypoallergenic Adhesive** The skin friendly hypoallergenic acrylic adhesive
 - Wound Dressing** To be used for wound protection and promotion of wound healing
 - Skin Barrier** An absorbable adhesive agent is used. It keeps the skin in the physiological condition.
 - Non-Adherent Pad** The absorptive part of the pad can not adhere to the wound surface.
 - Water-Proof / Drip-Proof** The water-proof and drip-proof structure prevents contamination from outside.

Blood Flow Care

More safety, easily, and comfortably —
Stagnation of the blood flow in the lower extremities during lying and rest is effectively prevented and the venous reflux is supported.

Corresponding to various needs and body shapes, the products line provides more appropriate venous return condition.

Medical Stocking
 (ANSILK-PRO J)



Calf Stocking



Stocking



Calf Stocking(Single Leg)
 Stocking(Single Leg)

100% Cotton Elastic Bandage
 (ELASCOT)





Compression & Gradation Care: Supporting venous return

Compression and Gradation care is a concept of care which ALCARE proposes. It can be achieved by compression and continuous changes of compression gradation.

Compression Design of ANSILK-PRO J

The three-dimensional structure fitting to the shape of leg and the compression design providing the optimum pressure according to Sigel theory* (38% increase in blood flow during rest) so it constantly gives adequate compression during lying.

Continuous Changes of Compression
 Venous Blood Flow




*Sigel et al, Arch Surg, Vol. 110 Feb. 1975

Quality of ANSILK-PRO J



Soft touch achieved without using latex

Color lines at the ankle and the knee can be a guide for your appropriate application. (Stocking type only)



Two-way delta stretch reduces the feeling of strong pressure on the thigh and assures comfortable fitting. (Stocking type only)



The "comfort heel" design increases the smooth feeling at the heel.



A monitor hole enables observation of blood flow conditions.

*Readjust the position (cover the toe) after observation. Leaving the toe uncovered may result in inadequate blood flow.

Compression & Gradation Care Products

■ Medical Stocking ■ Elastic Bandage

Stockings and bandages to support venous blood return before, during and after operation.

Medical Stocking



ANSILK-PRO J CALF STOCKINGS

Medical Stocking

- The calf stocking type causes no uncomfortable feeling and can be used by many people including those who usually wear no stockings.
- The material with soft texture assures comfortable fitting.
- The graduated compression is based on Sigel Theory.

Type	Product Code No.	Specification		Contents Per Box
		Girth of Ankle	Girth of Calf	
Both Legs	LL 17421	23~27cm	40~46cm	1 pair
	L 17422	21~25cm	36~42cm	
	M 17423	19~23cm	32~38cm	
	S 17424	17~21cm	28~34cm	
Single Leg	SS 17425	15~19cm	24~30cm	1 pce
	LL 17571	23~27cm	40~46cm	
	L 17572	21~25cm	36~42cm	
	M 17573	19~23cm	32~38cm	
	S 17574	17~21cm	28~34cm	
SS 17575	15~19cm	24~30cm		



ANSILK-PRO J STOCKINGS

Medical Stocking

- It rarely slips down from thigh because of the silicone garter.
- The graduated compression is based on Sigel Theory.
- The two-way delta stretch reduces feeling of strong pressure on the thigh.
- Color lines at the ankle and the knee can be a guide for your appropriate application.

Type	Product Code No.	Specification			Contents Per Box
		Girth of Ankle	Girth of Calf	Girth of Thigh	
Both Legs	LL 17431	23~27cm	40~46cm	55~64cm	1 pair
	L 17432	21~25cm	36~42cm	50~59cm	
	M 17433	19~23cm	32~38cm	47~56cm	
	S 17434	17~21cm	28~34cm	44~53cm	
Single Leg	SS 17435	15~19cm	24~30cm	40~49cm	1 pce
	LL 17581	23~27cm	40~46cm	55~64cm	
	L 17582	21~25cm	36~42cm	50~59cm	
	M 17583	19~23cm	32~38cm	47~56cm	
	S 17584	17~21cm	28~34cm	44~53cm	
SS 17585	15~19cm	24~30cm	40~49cm		

Elastic Bandage



ELASCOT

100% Cotton Elastic Bandage

- A 100% cotton elastic bandage has soft touch and excellent absorptivity.
- The original hard twist cotton yarn gives excellent elasticity to the bandage.
- It can be used as a substitute for medical stockings in the case that medical stocking is insufficient.

Type	Product Code No.	Specification	Contents Per Box
		Width×Length (stretched)	
2	11621	5.0cm×4.5m	6 rolls
3	11622	7.5cm×4.5m	6 rolls
4	11623	10.0cm×4.5m	6 rolls
7	11624	17.5cm×4.5m	6 rolls

Wound Care


The highly functional dressings as alternatives to gauze, the basic sanitary material, enable effective treatment, coverage and protection of wounds and promote patients' early recovery and efficiency of wound care.

Functions and Selection of Dressings



Feature Comparison of Dressings

● Optimum Features for After-Operative Suture Wounds
 ● Optimum Features for Wound Healing
 *Source: ALCARE's measurement data



Film Dressing
MULTI FIX
 Functions
 ● Air Permeability/ Water Proof
 ● Transparency


Easy

Wound Care Products

■ Hydrocolloid Dressing ■ Adhesive Island Dressing

ALCARE offers dressings from dry dressings to wet dressings, available for various uses such as wound protection and promotion of wound healing.

Hydrocolloid Dressing



KARAYAHESIVE

Hydrocolloid Wound Dressing

- Transparency enables follow-up observation of the wound.
- It can not slip even if it is applied to the joints or the parts with uneven surfaces because of appropriate flexibility and adherence.
- Appropriate absorbing exudate reduces maceration of normal skin.

EOG Sterilized Wound Dressing Water-Proof / Drip-Proof

Type	Product Code No.	Specification		Contents Per Box
		Size of Adhesive Part	Base Material	
C-3	15181	7.5cm×5.0cm	Urethane film	10 pcs
C-6	15182	15.0cm×5.0cm	Urethane film	10 pcs
C-9	15183	22.5cm×5.0cm	Urethane film	10 pcs
C-12	15184	30.0cm×5.0cm	Urethane film	10 pcs



KARAYAHESIVE

Hydrocolloid Wound Dressing

- The flesh color nonwoven is indistinctive on the face.
- The wet environment to promote wound healing is provided.
- Moderate adhesion causes minimum peeling irritation and assures easy exchange.

EOG Sterilized Wound Dressing

Type	Product Code No.	Specification		Contents Per Box
		Size of Adhesive Part	Base Material	
B-2	15191	5.0cm× 5.0cm	Polyester nonwoven	20 pcs
B-4	15192	10.0cm×10.0cm	Polyester nonwoven	10 pcs

Adhesive Island Dressing



MULTI FIX-PAD

Pad Attached Air Permeable / Waterproof Film Dressing

- The flexible backing film enables tight application without being wrinkled. Contamination from the outside and the effects of moisture are effectively blocked.
- Once the rayon pad absorbs the exudate, it curls to prevent its pad from adhering to the wound surface.
- The transparent adhesive part enables easy observation of the absorbed exudate.

EOG Sterilized Hypoallergenic Adhesive Non-Adherent Pad Water-Proof / Drip-Proof

Type	Product Code No.	Specification		Contents Per Box
		Size of Adhesive Part	Size of Absorptive Part	
2	16921	50mm× 65mm	27mm× 40mm	100 pcs
3	16922	80mm×100mm	40mm× 70mm	20 pcs
4	16923	100mm×125mm	55mm× 75mm	20 pcs
5	16924	80mm×150mm	40mm×110mm	20 pcs
6	16925	100mm×200mm	55mm×150mm	20 pcs
7	16926	100mm×250mm	55mm×200mm	20 pcs
8	16927	100mm×300mm	55mm×240mm	20 pcs
9	16928	100mm×370mm	55mm×310mm	20 pcs




SILKYPORE DRESSING

Adhesive Dressing for Wound Care

- A perforated mesh in combination with a non-woven pad allows uninterrupted breathing of the skin.
- A moderately soft dressing which can be applied to uneven areas of the skin.
- A gentle adhesive to prevent agitation to even delicate skins.

EOG Sterilized Hypoallergenic Adhesive Non-Adherent Pad

Type	Product Code No.	Specification		Contents Per Box
		Size of Adhesive Part	Size of Absorptive Part	
1	12008	4.0cm× 6.0cm	2.0cm× 3.5cm	100 pcs
1	15571	4.0cm× 6.0cm	2.0cm× 3.5cm	20 pcs
2	12001	5.0cm× 7.5cm	2.0cm× 4.0cm	50 pcs
2	15572	5.0cm× 7.5cm	2.0cm× 4.0cm	20 pcs
3	12002	6.0cm×10.0cm	3.0cm× 6.5cm	50 pcs
4	12003	10.0cm×13.0cm	4.0cm× 6.5cm	20 pcs
5	12004	10.0cm×16.0cm	5.0cm×10.0cm	20 pcs
6	12005	10.0cm×21.0cm	5.0cm×15.0cm	20 pcs
7	12006	10.0cm×26.0cm	5.0cm×20.0cm	20 pcs
8	12007	10.0cm×31.0cm	5.0cm×25.0cm	20 pcs



DERMAPORE DRESSING

Adhesive Dressing for Wound Care

- The air permeability and waterproofing property protects the wound from moisture and dirt.
- A moderately soft dressing which can be applied to uneven areas of the skin.
- A semi-transparent dressing allows inspection of the pad without the need for removal.

EOG Sterilized Hypoallergenic Adhesive Non-Adherent Pad Water-Proof / Drip-Proof

Type	Product Code No.	Specification		Contents Per Box
		Size of Adhesive Part	Size of Absorptive Part	
1	11991	4.0cm× 6.0cm	2.0cm× 3.5cm	100 pcs
1	15561	4.0cm× 6.0cm	2.0cm× 3.5cm	20 pcs
2	11992	5.0cm× 7.5cm	2.0cm× 4.0cm	50 pcs
2	15562	5.0cm× 7.5cm	2.0cm× 4.0cm	20 pcs
3	11993	6.0cm×10.0cm	3.0cm× 6.5cm	50 pcs
4	11994	10.0cm×13.0cm	4.0cm× 6.5cm	20 pcs
5	11995	10.0cm×16.0cm	5.0cm×10.0cm	20 pcs
6	11996	10.0cm×21.0cm	5.0cm×15.0cm	20 pcs
7	11997	10.0cm×26.0cm	5.0cm×20.0cm	20 pcs
8	11998	10.0cm×31.0cm	5.0cm×25.0cm	20 pcs

Wound Care Products

■ Absorbent Pad ■ Film Dressing ■ Fixation Material

Absorbent Pad



DERMA AID

Absorbent Pad for Wound

- 18 times the absorption power of gauze. Can be used on high exudates wounds.
- Non-adherent films, which rarely adhere to the wound surface, are applied to both sides of the pad.
- Each pad/pack is sterilized.

EOG Sterilized Non-Adherent Pad

Type	Product Code No.	Specification	Contents Per Box
1	11981	5.0cm× 5.0cm	100 pcs
2	11982	5.0cm× 7.5cm	100 pcs
4	11983	7.5cm×10.0cm	100 pcs
6	11984	10.0cm×10.0cm	100 pcs
8	11985	10.0cm×20.0cm	100 pcs

Film Dressing



MULTI FIX

Surgical Film with Air Permeability and Waterproofing Property

- The air permeability and waterproofing property do not permit penetration of water and external contamination.
- A transparent film enables easy follow-up observation of the wound.
- No.3 is large enough to cover KARAYAHESIVE C-3, while No.5 is for KARAYAHESIVE C-6.

EOG Sterilized Hypoallergenic Adhesive Water-Proof / Drip-Proof

Type	Product Code No.	Specification	Contents Per Box
1	16386	48mm× 48mm	100 pcs
2	16381	60mm× 70mm	100 pcs
3	16382	80mm×100mm	50 pcs
4	16383	100mm×125mm	50 pcs
5	16384	80mm×180mm	20 pcs
6	16385	100mm×200mm	20 pcs
S-2	16391	150mm×150mm	20 pcs

Fixation Material



DERMAPORE

Air Permeable Waterproof Elastic Adhesive Bandage

- The air permeability and waterproofing property do not permit penetration of water and external contamination.
- The extra thin bandage fits uneven sites because of its excellent elasticity in omni direction.
- A semi-transparent film enables easy follow-up observation of the wound.

Hypoallergenic Adhesive Water-Proof / Drip-Proof

Type	Product Code No.	Specification (W×L)	Contents Per Box
3	11951	2.5cm×10m	2 rolls
5	11952	5.0cm×10m	1 roll
7	11953	7.5cm×10m	1 roll
10	11954	10.0cm×10m	1 roll



SILKYPORE

Flexible Adhesive Gauze Bandage

- A perforated mesh in combination with a non-woven pad allows uninterrupted breathing of the skin.
- A moderately soft product which is easily applied to bending joints.
- Available in two colors to suit the application and the afflicted area.


Hypoallergenic Adhesive

Type	Product Code No.	Specification (W×L)	Contents Per Pack	
White	1	11921	1.25cm×10m	24 rolls
	3	11922	2.5cm×10m	12 rolls
	5	11923	5.0cm×10m	6 rolls
	7	11924	7.5cm×10m	4 rolls
	10	11925	10.0cm×10m	3 rolls
	15	11926	15.0cm×10m	2 rolls
Beige	30	11927	30.0cm×10m	1 roll
	3	11931	2.5cm×10m	12 rolls
	5	11932	5.0cm×10m	6 rolls
	7	11933	7.5cm×10m	4 rolls
10	11934	10.0cm×10m	3 rolls	

Drainage Care

The highly functional dressings as alternatives to gauze, the basic sanitary material, enable effective extracorporeal absorption and accumulation of drainage and promote patients' early recovery and efficiency of drainage care.

Function and Selection of the Products for Drainage Care




Drainage Care Products

■ Drainage Pouch ■ Drainage Pad ■ Absorbent Pad

ALCARE offers suitable products to support drainage care. Achieving better patients' QOL and efficient nursing treatment by reducing frequency of gauze exchange.

Drainage Pouch



SURGIDRAIN OPEN-TOP

Drainage after Operation

- By opening and closing the domed cap, it is possible to observe and treat the wound without removing the appliance.
- The domed cap can protect the wound from the external impact.
- A piece of drainage pouch can sufficiently manage 3 to 5 days drainage from the drain.

○ **γ-ray Sterilized** ○ **Skin Barrier** ○ **Water-Proof / Drip-Proof**

Type	Product Code No.	Specification		Contents Per Box
		Bag Size(L×W)	Max Available Size	
Post-Op-M	15671	205mm×130mm	50mm	2 pcs
S	15681	175mm×110mm	30mm	2 pcs
M	15682	205mm×130mm	50mm	2 pcs
L	15683	235mm×150mm	80mm×60mm	2 pcs
LL	15684	280mm×219mm	150mm×70mm	2 pcs



SURGIDRAIN ZIPPER

Drainage after Operation

- By opening and closing the zipper, it is possible to observe and treat the wound without removing the appliance.
- The flat form enables better fitting to body contours and reduces uncomfortable feeling during application.
- A piece of drainage pouch can sufficiently manage 3 to 5 days drainage from the drain.

○ **γ-ray Sterilized** ○ **Skin Barrier** ○ **Water-Proof / Drip-Proof**

Type	Product Code No.	Specification		Contents Per Box
		Bag Size(L×W)	Max Available Size	
S	15661	158mm×110mm	30mm	10 pcs
M	15662	188mm×130mm	50mm	10 pcs
L	15663	218mm×150mm	80mm×60mm	5 pcs



WELLCARE DRAIN


Fistula Drainage

- The thin flexible skin barrier for fistula has high pH buffering ability.
- The viscous drainage can be removed through the soft tap by milking the pouch.
- This product is available in three sizes to be applied to fistulae in different sizes.

○ **Skin Barrier** ○ **Water-Proof / Drip-Proof**

Type	Product Code No.	Specification		Contents Per Box
		Bag Size(L×W)	Max Available Size	
S	15831	160mm×95mm	35mm	5 pcs
M	15832	240mm×120mm	100mm×60mm	5 pcs
L	15833	280mm×160mm	120mm×100mm	5 pcs

Drainage Pad



SD PAD / SD PAD-WP

Drainage Pad

- The absorptivity is equivalent to about 40 gauze layers and its thickness is about half, which reduces uncomfortable feeling during application.
- A fibrous polymer absorber quickly absorbs drainage to reduce unpleasant feeling.
- Effective absorption reduces excessive moisture.

○ **γ-ray Sterilized** ○ **Water-Proof / Drip-Proof**

*Limited to WP(waterproof)type.

Type	Product Code No.	Specification		Contents Per Box
		Size	Weight Per Piece	
SD PAD 2	16751	150mm×150mm	Approx. 12g	50 pcs
SD PAD 3	16752	200mm×200mm	Approx. 21g	20 pcs
SD PAD-WP 3	16741	200mm×200mm	Approx. 20g	20 pcs
SD PAD-WP 4	16742	200mm×300mm	Approx. 30g	20 pcs

*Weight: Net weight of the pad excluding the sterilized bag.

*Do not apply this pad directly to the wound.

Absorbent Pad



DERMA AID

Absorbent Pad for Wound

- 18 times the absorption power of gauze. Can be used on high exudates wounds.
- Non-adhesive film on both sides of the dressing prevents damage to new cell growth when removed.
- Each pad/pack is sterilized.

○ **EOG Sterilized** ○ **Non-Adherent Pad**

Type	Product Code No.	Specification	Contents Per Box
1	11981	5.0cm×5.0cm	100 pcs
2	11982	5.0cm×7.5cm	100 pcs
4	11983	7.5cm×10.0cm	100 pcs
6	11984	10.0cm×10.0cm	100 pcs
8	11985	10.0cm×20.0cm	100 pcs

Drainage Care Products

■ Drainage Bag ■ Accessory ■ Skin Barrier for Post-PEG Treatment

Drainage Bag



DRAINAGE BAG

Bag for Drainage

- This bag can be connected to ALCARE's drainage pouches by simple handling and be used for accumulation of a large volume of drainage.
- A wide connector, which minimally narrows the tube diameter, assures smooth flow of drainage.

Type	Product Code No.	Specification		Contents Per Box
		Bag Size(L×W)	Capacity	
—	15841	320mm×230mm (tube length: 1m)	2000ml	5 pcs

Accessory



CONNECT CONE

Catheter Fixer

- Enable a secure fix of a catheter.
- Made from a flexible material which prevents the catheter from bending or twisting.
- By cutting the tip of the connector, it is possible to use catheters in the 6~40Fr range.

○ **EOG Sterilized**

Type	Product Code No.	Specification (Effective Catheter Size)	Contents Per Box
—	12061	6~40Fr	3 pcs



PROCARE SOFT WAFER STICK

Wafer Skin Barrier (Karaya Gum + Citrus Pectin)

- The stick-type skin barrier is enough flexible to be freely molded.
- This product assures the tight contact between the pouch and the uneven skin surface such as the peristoma skin area.

○ **Skin Barrier**

Type	Product Code No.	Specification (L×W)	Contents Per Box
Stick	15601	15mm×100mm	10 pcs



PROCARE POWDER

Powder Skin Barrier (Karaya Gum + Citrus Pectin)

- This is a mixed-type powder made of karaya gum and citrus pectin which are both excellent in their skin protecting quality.
- It fills the gap between the skin and the pouch to protect the skin.
- It is easy to apply, even to a small gap. The powder offers excellent water absorption.

○ **Skin Barrier**

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	12611	50g	1 bottle

Skin Barrier for Post-PEG Treatment



PEG CARE

Skin Barrier for Post-PEG Treatment

- This skin barrier is exclusively used in the treatment after percutaneous endoscopic gastrostomy (PEG).
- This thin skin barrier flexibly fits the skin without causing stress on the PEG kit.
- The skin barrier having the excellent pH buffering ability can effectively control exudates.

○ **Skin Barrier**

Type	Product Code No.	Specification (Outer Diameter)	Contents Per Box
—	16261	50mmφ	20 pcs

ALCARE's "Solution"

for all the people hoping for health

For the medical staff and the patients who seek for Best Medical Care, ALCARE takes a proactive stance on a broad range of issues in various clinical fields.

Pursuing of promotion of patient's early recovery and treatment efficiency

During / After Operation Management Solution



Hydrocolloid Wound Dressing <KARAYAHESIVE> Drainage after Operation <SURGIDRAIN ZIPPER> Pad Attached Air Permeable / Waterproof Film Dressing <MULTI FIX-PAD> Adhesive Dressing for Wound Care <SILKYPORE DRESSING>

Providing optimum materials for prevention and treatment

Wound Care Solution



Absorbent Pad for Wound <DERMA AID> Floatation Pad for Decubitus Prevention <BONMAT>

Pursuing a total support of blood flow care

Compression & Graduation Care Solution



Medical Stocking <ANSILK-II PANTY HOSE> 100% Cotton Elastic Bandage <ELASCOT> Pad for Compression <CG PAD> Medical Stocking <ANSILK-PRO J CALF STOCKINGS>

Realizing best care by providing a large variety of products

Bandage Care Solution



Spandex Elastic Bandage <HIGHSPAN> Elastic Net Bandage <PRESNET> Elastic Adhesive Bandage <ELATEX> Flexible Adhesive Gauze Bandage <SILKYPORE> Air Permeable Waterproof Elastic Adhesive Bandage <DERMAPORE>

Supporting improvement of ostomates' QOL


Stoma Care Solution



Open-Ended Bag for Colostomy / Ileostomy <LAPACK SUPER G> Open-Ended Bag for Colostomy / Ileostomy <POSPACK-K> One-Piece Ostomy System <YOU CARE D> Two-Piece Ostomy System <PROCARE 2 FA, PROCARE 2 C>

Improving patients' QOL by providing better continence care


Incontinence Care Solution



Male Incontinence Appliance <YUNIBON> Male Incontinence Appliance <LEG BAG>

Optimum products for the period from treatment to rehabilitation

Orthopedic Care Solution



Glass Fiber Casting Tape <CASTLIGHT-α> Plaster of Paris Bandage <PLASRUN> Water-Activated Glass Fiber Splint <LIGHT SPLINT-II> Servical Collar <POLYNECK HARD> Rib Belt <BUSTBAND-ACE> Lumbar Supporter <SACRO LIGHT-DX>

Protecting athletes from being injured

Sports Care Solution



<ANKLE GUARD EX> <LUMBAR GUARD> <KNEE GUARD> <REGUARD SKIN CARE>

Index

A	
ANSILK-PRO J CALF STOCKINGS.....	6
ANSILK-PRO J STOCKINGS.....	6
C	
CONNECT CONE.....	14
D	
DERMA AID.....	10,13
DERMAPORE.....	10
DERMAPORE DRESSING.....	9
DRAINAGE BAG.....	14
E	
ELASCOT.....	6
K	
KARAYAHESIVE.....	9
M	
MULTI FIX.....	10
MULTI FIX-PAD.....	9
P	
PEG CARE.....	14
PROCARE POWDER.....	14
PROCARE SOFT WAFER STICK.....	14
S	
SD PAD/SD PAD-WP.....	13
SILKYPORE.....	10
SILKYPORE DRESSING.....	9
SURGIDRAIN OPEN-TOP.....	13
SURGIDRAIN ZIPPER.....	13
W	
WELLCARE DRAIN.....	13

