

Bag Accessory

YOU CARE TRUNK BELT WITH FACE PLATE

Belt for One-Piece Bag




- Provides a securer application, preventing the flange from slipping or falling off with the body movement.
- It is recommendable for active athletes or for those who have difficulty in adhering the appliances to the skin.

Type	Product Code No.	Specification		Contents Per Pack
		Waist Size	Applicable Bag	
—	15081	60~110cm (diameter of hole 73mmφ)	YOU CARE TD YOU CARE D YOU CARE C YOU CARE U POSPACK-B POSPACK-LIGHT	1 set

POSPACK

Treatment for Immediately After Operation to Home-care



Simple & Easy

Simple and easy open-ended bag can be recommended for large volume stools.

After Operation & Home-care

Coverage from after operation to home-care with sterilization and skin barrier suitable for each condition.

ALCARE

ALCARE Co., Ltd.
INTERNATIONAL SALES DIVISION

19F, Arcacentral Bldg.
1-2-1 Kinshi, Sumida-ku
Tokyo 130-0013, Japan
TEL:81-3-5611-7819
FAX:81-3-5611-7827
<http://www.alcare.co.jp/e>

Printed in China
Sep 2007

ALCARE

POSPACK Series


Five types are available for various purposes, ranging from treatment immediately after operation to home-care.

Types and characteristics of POSPACK Series

Timing Exchange	Immediately after Operation		After Discharge (Home-care)	
	1-2 days	POSPACK-K STERILE EOG sterilized type is the optimum device to be used immediately after operation.	POSPACK-K POSPACK-K is a device with "KARAYA PLAST", a skin barrier which does not cause irritation.	POSPACK-LIGHT POSPACK-LIGHT is a device with an ultra-thin skin barrier that can be changed frequently. Pre-cut type is also available.
3-5 days			POSPACK-B POSPACK-B is a device that has enhanced adhesive strength and can be continuously used for long periods of time.	

POSPACK-B [Free-Cut]

Open-Ended Bag for Colostomy / Ileostomy



- The capacious open-ended bag with a large opening for removal is suitable for those who require greater capacity for accumulation of stool.

Type	Product Code No.	Specification		Contents Per Box
		Max Available Size	Ostomy Bag Size*	
—	15021	≤60mm (free-cut)	400mm×150mm	10 pcs

POSPACK-K STERILE [Free-Cut]

Open-Ended Bag for Colostomy / Ileostomy



- The EOG sterilized bag is the optimum appliance to be used immediately after operation.
- The open-ended bag with a large opening for removal facilitates enables easy peristomal treatment.
- Maximum acceptable stoma size up to 110 mm.

Type	Product Code No.	Specification		Contents Per Box
		Max Available Size	Ostomy Bag Size*	
70	15001	70mm	400mm×150mm	5 pcs
110	15002	110mm	400mm×170mm	5 pcs

EOG sterilized

POSPACK-LIGHT [Pre-Cut / Free-Cut]

Open-Ended Bag for Colostomy / Ileostomy



- A very thin gentle skin barrier.
- The capacious open-ended bag with a large opening for removal is suitable for those who require greater capacity for accumulation of stool.
- The free-cut type is also available.

Type	Product Code No.	Starter Hole Size	Specification		Contents Per Box
			Max Available Size	Ostomy Bag Size*	
20	15031	20mmφ	≤60mm (free-cut)	400mm×150mm	30 pcs
25	15032	25mmφ	19~24mm	400mm×150mm	30 pcs
30	15033	30mmφ	24~29mm	400mm×150mm	30 pcs
35	15034	35mmφ	29~34mm	400mm×150mm	30 pcs
40	15035	40mmφ	34~39mm	400mm×150mm	30 pcs
45	15036	45mmφ	39~44mm	400mm×150mm	30 pcs
50	15037	50mmφ	44~49mm	400mm×150mm	30 pcs

POSPACK-K [Free-Cut]

Open-Ended Bag for Colostomy / Ileostomy



- The open-ended bag with a large opening for removal facilitates enables easy peristomal treatment.
- Maximum acceptable stoma size up to 110 mm.

Type	Product Code No.	Specification		Contents Per Box
		Max Available Size	Ostomy Bag Size*	
70	15011	70mm	400mm×150mm	10 pcs
110	15012	110mm	400mm×170mm	10 pcs

POSPACK-SIMPLE [Free-Cut]

Open-Ended Bag for Colostomy / Ileostomy



- A very thin gentle skin barrier.
- The capacious open-ended bag with a large opening for removal is suitable for those who require greater capacity for accumulation of stool.
- The free-cut type is also available.
- This appliance is an economy type.

Type	Product Code No.	Specification		Contents Per Box
		Max Available Size	Ostomy Bag Size*	
—	15041	≤50mm (free-cut)	400mm×110mm	30 pcs

*Ostomy Bag Size : Length×Width

*Ostomy Bag Size : Length×Width