

ALCARE Ostomy Product Catalogue



History

Over 50 years with stoma care.

ALCARE remains committed to improving patients' quality of life.

1965 ● Launched "Lapack", the first adhesive ostomy appliance manufactured in Japan.



Lapack

1974 ● Launched "Karaya Disk", Japan's first skin barrier.

● Established "Happy Life Club" to support ostomates' independence.

1976 ● Launched "Happy Life Club News", an informative magazine to help ostomates to return to normal daily life.



First issue of
Happy Life Club News

1979 ● Launched "Pospack", a 1-piece ostomy appliance including a skin barrier.

1982 ● Launched "Biourine", a 1-piece urostomy appliance with a highly durable skin barrier.

1983 ● Held the first "Ostomates' Meeting" for ostomates, care workers, and the appliance manufacturer to exchange information with each other.



First Ostomates' Meeting

1985 ● Launched "Biomax", Japan's first 2-piece ostomy system.

1998 ● Launched "Youcare", a 1-piece ostomy appliance for short-term exchange.



Youcare

2001 ● Launched "Ilfine", a 1-piece ostomy appliance dedicated for ileostomy.



Ilfine

2008 ● Launched "Cellcare", a 2-piece ostomy system containing a skin care material.



Cellcare

INDEX

History	1
Digestive Ostomy Appliance	3
Youcare	4
Procare	8
Ilefine	11
Cellcare	12
Pospack	15
Lapack	15
Urinary Ostomy Appliance	17
Youcare	18
Procare	20
Cellcare	22
Accessory	24
Quick Guide for Accessory	34

Digestive Ostomy Appliance

Colostomy / Ileostomy

Description of icons

Faceplate (part to be attached to abdomen)

**Flat**

Flat and suitable for stoma with certain height.

**Safety Plate**

Soft plate included and good for both adhesion and followability.

**Convex Plate**

Convex plate included and suitable for stoma with not enough height.

Ostomy bag (the bag for collecting excrement)

**Open outlet type**

Open bag with outlet at bottom.

**Closed outlet type**

Closed bag without outlet.

**Cap outlet type**

Open bag with cap to close outlet for mushy and watery stool.

Free-cut**Free-cut**


Able to cut freely in accordance with a stoma size.

Pre-cut**Pre-cut**

No need to cut as a suitable hole is already made.

Youcare-TD 1-piece Ostomy System 2 days exchange

- Features**
- Elastic and durable thin skin barrier is used for faceplate.
 - Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool (Without-filter type is also available).
 - Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
Flesh	17441	20mmø	14~59mm(Free-cut)	275mmx150mm	10 Pieces
W/O filter	17451	20mmø	14~59mm(Free-cut)	275mmx150mm	10 Pieces
20	17461	20mmø	14~59mm(Free-cut)	275mmx150mm	10 Pieces
25	17462	25mmø	19~24mm	275mmx150mm	10 Pieces
30	17463	30mmø	24~29mm	275mmx150mm	10 Pieces
35	17464	35mmø	29~34mm	275mmx150mm	10 Pieces
40	17465	40mmø	34~39mm	275mmx150mm	10 Pieces
45	17466	45mmø	39~44mm	275mmx150mm	10 Pieces
50	17467	50mmø	44~49mm	275mmx150mm	10 Pieces

Youcare-TDc 1-piece Ostomy System 2 days exchange

- Features**
- Elastic and durable thin skin barrier is used for faceplate.
 - Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool.
 - Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
22	17471	22mmø	18~21mm	275mmx150mm	10 Pieces
25	17472	25mmø	21~24mm	275mmx150mm	10 Pieces
28	17473	28mmø	24~27mm	275mmx150mm	10 Pieces
32	17474	32mmø	27~31mm	275mmx150mm	10 Pieces
36	17475	36mmø	31~35mm	275mmx150mm	10 Pieces
40	17476	40mmø	35~39mm	275mmx150mm	10 Pieces
45	17477	45mmø	39~44mm	275mmx150mm	10 Pieces
50	17478	50mmø	44~49mm	275mmx150mm	10 Pieces

Youcare-D 1-piece Ostomy System 2 days exchange

- Features**
- Elastic and durable thin skin barrier is used for faceplate.
 - Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool (Without-filter type is also available).



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
Flesh	16761	20mmø	14~59mm(Free-cut)	275mmx150mm	30 Pieces
W/O filter	16771	20mmø	14~59mm(Free-cut)	275mmx150mm	30 Pieces
20	16781	20mmø	14~59mm(Free-cut)	275mmx150mm	30 Pieces
25	16782	25mmø	19~24mm	275mmx150mm	30 Pieces
30	16783	30mmø	24~29mm	275mmx150mm	30 Pieces
35	16784	35mmø	29~34mm	275mmx150mm	30 Pieces
40	16785	40mmø	34~39mm	275mmx150mm	30 Pieces
45	16786	45mmø	39~44mm	275mmx150mm	30 Pieces
50	16787	50mmø	44~49mm	275mmx150mm	30 Pieces

- Youcare Digestive
- Procare Digestive
- Ilefine Digestive
- Cellcare Digestive
- Others Digestive
- Youcare Urinary
- Procare Urinary
- Cellcare Urinary
- Accessory

<h2>Youcare-C</h2>	1-piece Ostomy System	1 day exchange
--------------------	-----------------------	----------------

Features

- Elastic ultra-thin skin barrier is used for faceplate.
- Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
Flesh	16791	20mmø	14~59mm(Free-cut)	200mm×150mm	30 Pieces
20	16801	20mmø	14~59mm(Free-cut)	200mm×150mm	30 Pieces
25	16802	25mmø	19~24mm	200mm×150mm	30 Pieces
30	16803	30mmø	24~29mm	200mm×150mm	30 Pieces
35	16804	35mmø	29~34mm	200mm×150mm	30 Pieces
40	16805	40mmø	34~39mm	200mm×150mm	30 Pieces
45	16806	45mmø	39~44mm	200mm×150mm	30 Pieces
50	16807	50mmø	44~49mm	200mm×150mm	30 Pieces


Youcare 2-F

2-piece Ostomy System

3-5 days exchange

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Thinner outer edge of faceplate assures close contact with body and reduces uncomfortableness while attached.



Free-cut

Pre-cut

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M22	16591	22mmø	16~34mm(Free-cut)	97mm×105mm	M	5 Pieces
M25	16592	25mmø	21~24mm	97mm×105mm	M	5 Pieces
M28	16593	28mmø	24~27mm	97mm×105mm	M	5 Pieces
L32	16594	32mmø	27~31mm	108mm×117mm	L	5 Pieces
L36	16595	36mmø	31~35mm	108mm×117mm	L	5 Pieces
L40	16596	40mmø	35~39mm	108mm×117mm	L	5 Pieces
L45	16597	10mmø	4~44mm(Free-cut)	108mm×117mm	L	5 Pieces
65	16598	10mmø	4~64mm(Free-cut)	129mm×140mm	65	5 Pieces

Youcare 2-Fc

2-piece Ostomy System

3-5 days exchange

Features

- Tape around faceplate prevents skin barrier from dissolving while bathing.



Pre-cut

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M22	16601	22mmø	18~21mm	110mmø	M	5 Pieces
M25	16602	25mmø	21~24mm	110mmø	M	5 Pieces
M28	16603	28mmø	24~27mm	110mmø	M	5 Pieces
L32	16604	32mmø	27~31mm	120mmø	L	5 Pieces
L36	16605	36mmø	31~35mm	120mmø	L	5 Pieces
L40	16606	40mmø	35~39mm	120mmø	L	5 Pieces

Youcare 2-TDf

2-piece Ostomy System

Features

- Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool.
- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.
- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification		Contents Per Box
		Bag Size(L×W)		
Transparent	M	17761	300mm×150mm	10 Pieces
	L	17762	300mm×150mm	10 Pieces
	65	17763	300mm×150mm	10 Pieces
Flesh	M	17771	300mm×150mm	10 Pieces
	L	17772	300mm×150mm	10 Pieces

Youcare Digestive

Youcare 2-TD

2-piece Ostomy System

Features

- Stoma bag without a deodorizing filter.
- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.
- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification		Contents Per Box
		Bag Size(LxW)		
Transparent	M	17781	300mm×150mm	10 Pieces
	L	17782	300mm×150mm	10 Pieces
	65	17783	300mm×150mm	10 Pieces
Flesh	M	17791	300mm×150mm	10 Pieces
	L	17792	300mm×150mm	10 Pieces

Procure Digestive

Ilefine Digestive

Cellcare Digestive

Youcare 2-Df

2-piece Ostomy System

Features

- Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool.
- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification		Contents Per Box
		Bag Size(LxW)		
Transparent	M	16611	300mm×150mm	10 Pieces
	L	16612	300mm×150mm	10 Pieces
	65	16613	300mm×150mm	10 Pieces
Flesh	M	17741	300mm×150mm	10 Pieces
	L	17742	300mm×150mm	10 Pieces

Others Digestive

Youcare Urinary

Procure Urinary

Cellcare Urinary

Youcare 2-D

2-piece Ostomy System

Features

- Stoma bag without a deodorizing filter.
- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification		Contents Per Box
		Bag Size(LxW)		
Transparent	M	16621	300mm×150mm	10 Pieces
	L	16622	300mm×150mm	10 Pieces
	65	16623	300mm×150mm	10 Pieces
Flesh	M	17751	300mm×150mm	10 Pieces
	L	17752	300mm×150mm	10 Pieces

Accessory

Procare 1-Free-Cut D

1-piece Ostomy System

3 to 5 days exchange

Features

- Comfortable non-woven fabric is attached to both sides of stoma bag.
- Can be attached while checking position of stoma from slit of outer non-woven fabric.
- Tape around faceplate prevents skin barrier from dissolving while bathing.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
50	12331	20mmø	14~49mm(Free-cut)	300mm×150mm	10 Pieces



Free-cut

Procare 1-Pre-Cut Dc

1-piece Ostomy System

3 to 5 days exchange

Features

- Tape around faceplate prevents skin barrier from dissolving while bathing.
- Trunk belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
22	12366	22mmø	18~21mm	300mm×150mm	5 Pieces
25	12361	25mmø	21~24mm	300mm×150mm	5 Pieces
28	12362	28mmø	24~27mm	300mm×150mm	5 Pieces
32	12363	32mmø	27~31mm	300mm×150mm	5 Pieces
36	12364	36mmø	31~35mm	300mm×150mm	5 Pieces
40	12365	40mmø	35~39mm	300mm×150mm	5 Pieces



Pre-cut

Procare 1-Free-Cut C

1-piece Ostomy System

1 day exchange

Features

- Comfortable non-woven fabric is attached to both sides of stoma bag.
- Can be attached while checking position of stoma from slit of outer non-woven fabric.
- Tape around faceplate prevents skin barrier from dissolving while bathing.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
50	12341	20mmø	14~49mm(Free-cut)	210mm×150mm	10 Pieces



Free-cut

Youcare Digestive

Procare Digestive

Ilefine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procare Urinary

Cellcare Urinary

Accessory

Pediatric Procare 1-Post Op

1-piece Ostomy System

1 day exchange

Features

- Open-ended bag with wide outlet for removal makes peristomal care easier.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
—	14211	—	~34mm(Free-cut)	230mm×110mm	10 Pieces



Free-cut

Pediatric Procare 1-D

1-piece Ostomy System

2 to 3 days exchange

Features

- Oval faceplate assures sufficient adhesive area while not interfering with body movements.
- Pediatric Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
—	14221	—	24×34mm(Free-cut)	230mm×130mm	10 Pieces



Free-cut

Procare 2-F

2-piece Ostomy System

3 to 5 days exchange

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Available stoma size: Up to a maximum diameter of 69mm.



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
30	12091	10mmø	4~29mm(Free-cut)	100mm×100mm	30	5 Pieces
40	12092	10mmø	29~39mm(Free-cut)	100mm×100mm	40	5 Pieces
50	12093	10mmø	39~49mm(Free-cut)	110mm×110mm	50	5 Pieces
70	12094	10mmø	49~69mm(Free-cut)	130mm×130mm	70	5 Pieces



Free-cut


Procare 2-FA

2-piece Ostomy System

3 to 5 days exchange

Features

- Tape around faceplate prevents skin barrier from dissolving while bathing.



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
30	12101	10mmø	4~29mm(Free-cut)	100mmø	30	5 Pieces
40	12102	10mmø	29~39mm(Free-cut)	110mmø	40	5 Pieces
50	12103	10mmø	39~49mm(Free-cut)	120mmø	50	5 Pieces



Free-cut


Procare 2-Fc

2-piece Ostomy System

3 to 5 days exchange

Features

- Tape around faceplate prevents skin barrier from dissolving while bathing.
- Trunk belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
22-40	12126	22mmø	18~21mm	110mmø	40	5 Pieces
25-40	12121	25mmø	21~24mm	110mmø	40	5 Pieces
28-50	12122	28mmø	24~27mm	120mmø	50	5 Pieces
32-50	12123	32mmø	27~31mm	120mmø	50	5 Pieces
36-50	12124	36mmø	31~35mm	120mmø	50	5 Pieces
40-50	12125	40mmø	35~39mm	120mmø	50	5 Pieces



Pre-cut

Procare 2-D

2-piece Ostomy System

Features

- Transparent bag allows you to observe stool condition and flesh bag is less noticeable.



Type	Product Code No.	Specification		Contents Per Box
		Bag Size(LxW)		
Transparent	30	12131	300mmx150mm	10 Pieces
	40	12132	300mmx150mm	10 Pieces
	50	12133	300mmx150mm	10 Pieces
	70	12134	310mmx165mm	10 Pieces
Flesh	30	13641	300mmx150mm	10 Pieces
	40	13642	300mmx150mm	10 Pieces
	50	13643	300mmx150mm	10 Pieces



Youcare Digestive

Procure Digestive

Ilefine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procure Urinary

Cellcare Urinary

Accessory

Ilefine Special model for Ileostomy For digestive stoma

Ilefine

**Ostomy Appliance
for Ileostomy**


2 to 3 days exchange

Features

- Faceplate has skin barrier with excellent pH buffering action which can protect skin from strong alkaline stool.
 - Faceplate with Safety Plate which prevents watery stool from sneaking in.
- *D-70 and Dcap flat: Flat faceplate without Safety Plate.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
D-Cap 30	16991	15mmø	14~29mm(Free-cut)	346mm×148mm	10 Pieces
D-Cap 40	16992	25mmø	24~39mm(Free-cut)	346mm×148mm	10 Pieces
D-Cap 50	16993	35mmø	34~49mm(Free-cut)	346mm×148mm	10 Pieces
D-Cap Flat 40	18171	15mmø	14~39mm(Free-cut)	346mm×148mm	10 Pieces
D-Cap Flat 50	18172	15mmø	14~49mm(Free-cut)	346mm×148mm	10 Pieces
D-Cap Flat 60	18173	15mmø	14~59mm(Free-cut)	346mm×148mm	10 Pieces
TD-30	19581	15mmø	14~29mm(Free-cut)	300mm×148mm	10 Pieces
TD-40	19582	25mmø	24~39mm(Free-cut)	300mm×148mm	10 Pieces
TD-50	19583	35mmø	34~49mm(Free-cut)	300mm×148mm	10 Pieces
D-30	16981	15mmø	14~29mm(Free-cut)	300mm×148mm	10 Pieces
D-40	16982	25mmø	24~39mm(Free-cut)	300mm×148mm	10 Pieces
D-50	16983	35mmø	34~49mm(Free-cut)	300mm×148mm	10 Pieces
D-70	16984	20mmø	14~69mm(Free-cut)	300mm×148mm	5 Pieces



Youcare Digestive

Procare Digestive

Ilefine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procare Urinary

Cellcare Urinary

Accessory

Cellcare 1-TD

1-piece Ostomy System

2 to 5 days exchange

Features

- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.
- “Air permeability recovery filter”, which helps the filter to be re-functioned when clogged, prevents the leakage of odor and stool.
- Reduces uncomfortableness while attached as faceplate and bag are put together.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
Flesh	19421	20mmø	14~49mm(Free-cut)	300mmx150mm	10 Pieces
20 Free	18701	20mmø	14~49mm(Free-cut)	300mmx150mm	10 Pieces
25	18702	25mmø	19~24mm	300mmx150mm	10 Pieces
30	18703	30mmø	24~29mm	300mmx150mm	10 Pieces
35	18704	35mmø	29~34mm	300mmx150mm	10 Pieces
40	18705	40mmø	34~39mm	300mmx150mm	10 Pieces
70 Free	18706	10mmø	4~69mm(Free-cut)	300mmx150mm	5 Pieces

Cellcare 1-TDs

1-piece Ostomy System

2 to 5 days exchange

Features

- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.
- “Air permeability recovery filter”, which helps the filter to be re-functioned when clogged, prevents the leakage of odor and stool.
- Reduces uncomfortableness while attached as faceplate and bag are put together.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
30 Free	20151	15mmø	14~29mm(Free-cut)	300mmx150mm	5 Pieces
40 Free	20152	25mmø	24~39mm(Free-cut)	300mmx150mm	5 Pieces
50 Free	20153	35mmø	34~49mm(Free-cut)	300mmx150mm	5 Pieces

Cellcare 1-TDc

1-piece Ostomy System

2 to 5 days exchange

Features

- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.
- “Air permeability recovery filter”, which helps the filter to be re-functioned when clogged, prevents the leakage of odor and stool.
- Reduces uncomfortableness while attached as faceplate and bag are put together.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
22	18711	22mmø	18~21mm	300mmx150mm	5 Pieces
25	18712	25mmø	21~24mm	300mmx150mm	5 Pieces
28	18713	28mmø	24~27mm	300mmx150mm	5 Pieces
32	18714	32mmø	27~31mm	300mmx150mm	5 Pieces
36	18715	36mmø	31~35mm	300mmx150mm	5 Pieces
40	18716	40mmø	35~39mm	300mmx150mm	5 Pieces
45	18717	45mmø	39~44mm	300mmx150mm	5 Pieces
50	18718	50mmø	44~49mm	300mmx150mm	5 Pieces
45 Free	20141	15mmø	29~44mm(Free-cut)	300mmx150mm	5 Pieces
55 Free	20142	25mmø	39~54mm(Free-cut)	300mmx150mm	5 Pieces

Youcare Digestive

Procare Digestive

Ilefine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procare Urinary

Cellcare Urinary

Accessory

Cellcare 1-D Cap

1-piece Ostomy System

2 to 5 days exchange

Features

- Soft type outlet enables discharge of watery or muddy stool by milking.
- "Air permeability recovery filter", which helps the filter to be re-functioned when clogged, prevents the leakage of odor and stool.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
Flesh	19431	20mmø	14~49mm(Free-cut)	346mmx148mm	10 Pieces
20 Free	18921	20mmø	14~49mm(Free-cut)	346mmx148mm	10 Pieces
25	18922	25mmø	19~24mm	346mmx148mm	10 Pieces
30	18923	30mmø	24~29mm	346mmx148mm	10 Pieces
35	18924	35mmø	29~34mm	346mmx148mm	10 Pieces
40	18925	40mmø	34~39mm	346mmx148mm	10 Pieces
70 Free	18926	10mmø	4~69mm(Free-cut)	346mmx148mm	10 Pieces

Cellcare 1-Dc Cap

1-piece Ostomy System

2 to 5 days exchange

Features

- Soft type outlet enables discharge of watery or muddy stool by milking.
- "Air permeability recovery filter", which helps the filter to be re-functioned when clogged, prevents the leakage of odor and stool.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
22	18931	22mmø	18~21mm	346mmx148mm	5 Pieces
25	18932	25mmø	21~24mm	346mmx148mm	5 Pieces
28	18933	28mmø	24~37mm	346mmx148mm	5 Pieces
32	18934	32mmø	27~31mm	346mmx148mm	5 Pieces
36	18935	36mmø	31~35mm	346mmx148mm	5 Pieces
40	18936	40mmø	35~39mm	346mmx148mm	5 Pieces
45 Free	18937	15mmø	29~44mm(Free-cut)	346mmx148mm	5 Pieces
55 Free	18938	25mmø	39~54mm(Free-cut)	346mmx148mm	5 Pieces

Cellcare 2-F

2-piece Ostomy System

2 to 5 days exchange

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Soft and thin skin barrier fits skin well.



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M Free	18251	10mmø	4~39mm(Free-cut)	88mmx88mm	M	5 Pieces
M16	18661	16mmø	12~15mm	88mmx88mm	M	5 Pieces
M19	18662	19mmø	15~18mm	88mmx88mm	M	5 Pieces
M22	18252	22mmø	18~21mm	88mmx88mm	M	5 Pieces
M25	18253	25mmø	21~24mm	88mmx88mm	M	5 Pieces
M28	18254	28mmø	24~27mm	88mmx88mm	M	5 Pieces
L Free	18255	10mmø	4~49mm(Free-cut)	98mmx98mm	L	5 Pieces
L32	18256	32mmø	27~31mm	98mmx98mm	L	5 Pieces
L36	18257	36mmø	31~35mm	98mmx98mm	L	5 Pieces
L40	18258	40mmø	35~39mm	98mmx98mm	L	5 Pieces
LL Free	18259	10mmø	4~69mm(Free-cut)	118mmx118mm	LL	5 Pieces

Cellcare 2-Fc

2-piece Ostomy System

2 to 5 days exchange

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Soft faceplate fits skin well.
- Trunk Belt II can improve close contact between faceplate and skin.



Pre-cut

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M13	18387	13mmø	9~12mm	100mmø	M	5 Pieces
M16	18388	16mmø	12~15mm	100mmø	M	5 Pieces
M19	18389	19mmø	15~18mm	100mmø	M	5 Pieces
M22	18381	22mmø	18~21mm	100mmø	M	5 Pieces
M25	18382	25mmø	21~24mm	100mmø	M	5 Pieces
M28	18383	28mmø	24~27mm	100mmø	M	5 Pieces
L32	18384	32mmø	27~31mm	110mmø	L	5 Pieces
L36	18385	36mmø	31~35mm	110mmø	L	5 Pieces
L40	18386	40mmø	35~39mm	110mmø	L	5 Pieces

Cellcare 2-TDf

2-piece Ostomy System

Features

- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.
- "Air permeability recovery filter", which helps the filter to be re-functioned when clogged, prevents the leakage of odor and stool.
- Hygienic and clean because stool adhering around joint ring of faceplate can be removed when stoma bag is exchanged.



Type	Product Code No.	Specification	Contents Per Box	
		Bag Size(LxW)		
Transparent	M	18261	300mmx150mm	10 Pieces
	L	18262	300mmx150mm	10 Pieces
	LL	18263	300mmx150mm	10 Pieces
Flesh	M	18671	300mmx150mm	10 Pieces
	L	18672	300mmx150mm	10 Pieces

Cellcare 2-D Cap

2-piece Ostomy System

Features

- Soft type outlet enables discharge of watery or muddy stool by milking.
- "Air permeability recovery filter", which helps the filter to be re-functioned when clogged, prevents the leakage of odor and stool.



Type	Product Code No.	Specification	Contents Per Box	
		Bag Size(LxW)		
Transparent	M	18911	346mmx148mm	10 Pieces
	L	18912	346mmx148mm	10 Pieces
	LL	18913	346mmx148mm	10 Pieces
Flesh	M	19441	346mmx148mm	10 Pieces
	L	19442	346mmx148mm	10 Pieces

Youcare Digestive

Procare Digestive

Ilefine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procare Urinary

Cellcare Urinary

Accessory

Youcare Digestive

Procare Digestive

Ilefine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procare Urinary

Cellcare Urinary

Accessory

Pospack-B

Open-Ended Bag for Colostomy/Ileostomy

3 to 5 days exchange

Features

- Recommended for people who have a large amount of stool thanks to open-ended bag with large capacity and wide outlet for removal.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
—	15021	20mmø	14~59mm(Free-cut)	400mmx150mm	10 Pieces



Free-cut

Pospack-Simple

Open-Ended Bag for Colostomy/Ileostomy

1 day exchange

Features

- Elastic ultra-thin skin barrier is used for faceplate.
- Recommended for people who have a large amount of stool thanks to open-ended bag with large capacity and wide outlet for removal.
- Very economical and cost effective.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
—	15041	20mmø	14~49mm(Free-cut)	400mmx110mm	30 Pieces



Free-cut

Lapack-Super G

Open-Ended Bag for Colostomy/Ileostomy

Features

- Open-ended bag allows you to dispose stool easily through wide outlet.
- Recommended for people who have a large amount of stool thanks to open-ended bag with large capacity and wide outlet for removal.
- Very economical and cost effective.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
80	12501	—	~39mm(Free-cut)	300mmx90mm	30 Pieces/Pack
100	12502	—	~49mm(Free-cut)	400mmx110mm	30 Pieces/Pack
125	12503	—	~74mm(Free-cut)	400mmx135mm	30 Pieces/Pack



Free-cut

Lapack-Super C

Closed-Ended Bag for Colostomy

Features

- Compact and light weight.
- Easy for application and observation thanks to transparent bag.
- Very economical and cost effective.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
100	12511	30mmø	24~49mm(Free-cut)	200mmx140mm	30 Pieces/Pack
125	12512	30mmø	24~74mm(Free-cut)	250mmx150mm	30 Pieces/Pack



Youcare Digestive

Procare Digestive

Ilefine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procare Urinary

Cellcare Urinary

Accessory

Urinary Ostomy Appliance

Urostomy

Description of icons

Faceplate (part to be attached to abdomen)



Flat

Flat and suitable for stoma with certain height.



Convex Plate

Convex plate included and suitable for stoma with not enough height.

Ostomy bag (the bag for collecting excrement)



Urine outlet type (Double lock)

Urostomy bag with safety double lock to close outlet.



Urine outlet type (Single lock)

Urostomy bag with easy single lock to close outlet.

Free-cut

Free-cut

Able to cut freely in accordance with a stoma size.

Pre-cut

Pre-cut

No need to cut as a suitable hole is already made.

Youcare-U

1-piece Ostomy System

2 to 3 days exchange

Features

- Elastic and durable thin skin barrier, which is specially developed for the urostomy, is used for faceplate.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
10	15361	10mmø	4~59mm(Free-cut)	260mm×150mm	10 Pieces
20	15362	20mmø	14~19mm	260mm×150mm	10 Pieces
25	15363	25mmø	19~24mm	260mm×150mm	10 Pieces
30	15364	30mmø	24~29mm	260mm×150mm	10 Pieces
35	15365	35mmø	29~34mm	260mm×150mm	10 Pieces
40	15366	40mmø	34~39mm	260mm×150mm	10 Pieces


Youcare-Uc

1-piece Ostomy System

2 to 3 days exchange

Features

- Elastic and durable thin skin barrier, which is specially developed for the urostomy, is used for faceplate.
- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
13	16401	13mmø	9~12mm	260mm×150mm	10 Pieces
16	16402	16mmø	12~15mm	260mm×150mm	10 Pieces
19	16403	19mmø	15~18mm	260mm×150mm	10 Pieces
22	16404	22mmø	18~21mm	260mm×150mm	10 Pieces
25	16405	25mmø	21~24mm	260mm×150mm	10 Pieces
28	16406	28mmø	24~27mm	260mm×150mm	10 Pieces
32	16407	32mmø	27~31mm	260mm×150mm	10 Pieces

Youcare 2-Fu

2-piece Ostomy System

4 to 7 days exchange

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Thinner outer edge of faceplate assures close contact with body and reduces uncomfortableness while attached.



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M13	16641	13mmø	7~34mm(Free-cut)	97mm×105mm	M	5 Pieces
M16	16642	16mmø	12~15mm	97mm×105mm	M	5 Pieces
M19	16643	19mmø	15~18mm	97mm×105mm	M	5 Pieces
M22	16644	22mmø	18~21mm	97mm×105mm	M	5 Pieces
M25	16645	25mmø	21~24mm	97mm×105mm	M	5 Pieces
M28	16646	28mmø	24~27mm	97mm×105mm	M	5 Pieces
L32	16647	32mmø	27~44mm(Free-cut)	108mm×117mm	L	5 Pieces

Youcare
Digestive

Procare
Digestive

Ilefine
Digestive

Cellcare
Digestive

Others
Digestive

Youcare
Urinary

Procare
Urinary

Cellcare
Urinary

Accessory

Youcare 2-Fuc

2-piece Ostomy System

4 to 7 days exchange

Features

- Tape around faceplate prevents skin barrier from dissolving while bathing.



Pre-cut

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M13	16651	13mmø	9~12mm	110mmø	M	5 Pieces
M16	16652	16mmø	12~15mm	110mmø	M	5 Pieces
M19	16653	19mmø	15~18mm	110mmø	M	5 Pieces
M22	16654	22mmø	18~21mm	110mmø	M	5 Pieces
M25	16655	25mmø	21~24mm	110mmø	M	5 Pieces
M28	16656	28mmø	24~27mm	110mmø	M	5 Pieces
L32	16657	32mmø	27~31mm	120mmø	L	5 Pieces

Youcare 2-U

2-piece Ostomy System

Features

- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification	Contents Per Box
		Bag Size(L×W)	
M	16661	260mm×150mm	10 Pieces
L	16662	260mm×150mm	10 Pieces

Procare 1-Free-Cut U

1-piece Ostomy System

4 to 7 days exchange

Features

- Comfortable non-woven fabric is attached to both sides of stoma bag.
- Can be attached while checking position of stoma from slit of outer non-woven fabric.
- Tape around faceplate prevents skin barrier from dissolving while bathing.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
50	12351	10mmø	4~49mm(Free-cut)	260mm×150mm	10 Pieces



Free-cut

Procare 1-Pre-Cut Uc

1-piece Ostomy System

4 to 7 days exchange

Features

- Tape around faceplate prevents skin barrier from dissolving while bathing.
- Trunk belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
16	12371	16mmø	12~15mm	260mm×150mm	5 Pieces
19	12372	19mmø	15~18mm	260mm×150mm	5 Pieces
22	12373	22mmø	18~21mm	260mm×150mm	5 Pieces
25	12374	25mmø	21~24mm	260mm×150mm	5 Pieces
28	12375	28mmø	24~27mm	260mm×150mm	5 Pieces



Pre-cut

Pediatric Procure 1-U

1-piece Ostomy System

2 to 3 days exchange

Features

- Oval faceplate assures sufficient adhesive area while not interfering with body movements.
- Pediatric Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
—	14231	—	24mm×34mm(Free-cut)	205mm×130mm	10 Pieces



Free-cut

Youcare
Digestive

Procare
Digestive

Ilefine
Digestive

Cellcare
Digestive

Others
Digestive

Youcare
Urinary

Procare
Urinary

Cellcare
Urinary

Accessory


Procare 2-FA Uro

2-piece Ostomy System

4 to 7 days exchange

Features

- Tape around faceplate prevents skin barrier from dissolving while bathing.



Free-cut

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
30	13381	10mmø	4~29mm(Free-cut)	100mmø	30	5 Pieces
40	13382	10mmø	29~39mm(Free-cut)	110mmø	40	5 Pieces
50	13383	10mmø	39~49mm(Free-cut)	120mmø	50	5 Pieces

Procare 2-Fc Uro

2-piece Ostomy System

4 to 7 days exchange

Features

- Tape around faceplate prevents skin barrier from dissolving while bathing.
- Trunk belt can improve close contact between faceplate and skin.



Pre-cut

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
16-40	13391	16mmø	12~15mm	110mmø	40	5 Pieces
19-40	13392	19mmø	15~18mm	110mmø	40	5 Pieces
22-40	13393	22mmø	18~21mm	110mmø	40	5 Pieces
25-40	13394	25mmø	21~24mm	110mmø	40	5 Pieces
28-50	13395	28mmø	24~27mm	120mmø	50	5 Pieces

Procare 2-U

2-piece Ostomy System

Features

- Transparent bag allows you to observe urine condition and flesh bag is less noticeable.



Type	Product Code No.	Specification		Contents Per Box
		Bag Size(LxW)		
Transparent	30	12151	260mm×150mm	10 Pieces
	40	12152	260mm×150mm	10 Pieces
	50	12153	260mm×150mm	10 Pieces
Flesh	30	13361	260mm×150mm	10 Pieces
	40	13362	260mm×150mm	10 Pieces
	50	13363	260mm×150mm	10 Pieces

Cellcare 1-U

1-piece Ostomy System

2 to 5 days exchange

Features

- Reduces uncomfortableness while attached as faceplate and bag are put together.
- Double lock outlet reduces concern for leakage.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
10 Free	18721	10mmø	4~49mm(Free-cut)	260mmx150mm	10 Pieces
20	18722	20mmø	14~19mm	260mmx150mm	10 Pieces
25	18723	25mmø	19~24mm	260mmx150mm	10 Pieces
30	18724	30mmø	24~29mm	260mmx150mm	10 Pieces
35	18725	35mmø	29~34mm	260mmx150mm	10 Pieces
40	18726	40mmø	34~39mm	260mmx150mm	10 Pieces

Cellcare 1-Uc

1-piece Ostomy System

2 to 5 days exchange

Features

- Convex faceplate helps to make sufficient height for stoma and prevents leakage of stool.
- Double lock outlet reduces concern for leakage.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
13	18731	13mmø	9~12mm	260mmx150mm	5 Pieces
16	18732	16mmø	12~15mm	260mmx150mm	5 Pieces
19	18733	19mmø	15~18mm	260mmx150mm	5 Pieces
22	18734	22mmø	18~21mm	260mmx150mm	5 Pieces
25	18735	25mmø	21~24mm	260mmx150mm	5 Pieces
28	18736	28mmø	24~27mm	260mmx150mm	5 Pieces
32	18737	32mmø	27~31mm	260mmx150mm	5 Pieces

Cellcare 2-F

2-piece Ostomy System

2 to 5 days exchange

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Soft and thin skin barrier fits skin well.



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M Free	18251	10mmø	4~39mm(Free-cut)	88mmx88mm	M	5 Pieces
M16	18661	16mmø	12~15mm	88mmx88mm	M	5 Pieces
M19	18662	19mmø	15~18mm	88mmx88mm	M	5 Pieces
M22	18252	22mmø	18~21mm	88mmx88mm	M	5 Pieces
M25	18253	25mmø	21~24mm	88mmx88mm	M	5 Pieces
M28	18254	28mmø	24~27mm	88mmx88mm	M	5 Pieces
L Free	18255	10mmø	4~49mm(Free-cut)	98mmx98mm	L	5 Pieces
L32	18256	32mmø	27~31mm	98mmx98mm	L	5 Pieces
L36	18257	36mmø	31~35mm	98mmx98mm	L	5 Pieces
L40	18258	40mmø	35~39mm	98mmx98mm	L	5 Pieces

Youcare Digestive

Procare Digestive

Ilfine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procare Urinary

Cellcare Urinary

Accessory

Cellcare 2-Fc

2-piece Ostomy System

2 to 5 days exchange

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Soft faceplate fits skin well.
- Trunk Belt II can improve close contact between faceplate and skin.



Pre-cut

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M13	18387	13mmø	9~12mm	100mmø	M	5 Pieces
M16	18388	16mmø	12~15mm	100mmø	M	5 Pieces
M19	18389	19mmø	15~18mm	100mmø	M	5 Pieces
M22	18381	22mmø	18~21mm	100mmø	M	5 Pieces
M25	18382	25mmø	21~24mm	100mmø	M	5 Pieces
M28	18383	28mmø	24~27mm	100mmø	M	5 Pieces
L32	18384	32mmø	27~31mm	110mmø	L	5 Pieces
L36	18385	36mmø	31~35mm	110mmø	L	5 Pieces
L40	18386	40mmø	35~39mm	110mmø	L	5 Pieces

Cellcare 2-U

2-piece Ostomy System

Features

- Double lock outlet reduces concern for leakage.



Type	Product Code No.	Specification	Contents Per Box
		Bag Size(L×W)	
M	18681	260mm×150mm	10 Pieces
L	18682	260mm×150mm	10 Pieces

Accessory

Procare Remover

Remover

Features

- Reduces irritation upon peeling ostomy appliance.
- Can also be used for immediate removal when leakage etc..
- Easy to use and economical thanks to individual packing.



Type	Product Code No.	Specification(Liquid Amount)	Contents Per Box
—	13871	2mL	50 Pieces

Remois Cleanse

Skin Moisture
Cleansing Cream

Features

- Natural oil ingredient detaches dirt from skin and wipe it off for clean skin.
- Fresh feeling sense of use and ostomy appliance can be applied right after wiping off.
- Moisturizing ingredient keeps skin moist and smooth and protects it from drying.



Type	Product Code No.	Specification(Contents Amount)	Contents Per Box
Regular	17701	180g	1 Piece
Push Bottle	17703	500g	1 Piece
Handy	17702	5g/Piece	10 Pieces

Procure Powder

Powder Skin Barrier

Features

- Protects skin from excrement by filling gaps between stoma and faceplate.
- Easy to apply even to small gaps and has excellent water absorption.



Type	Product Code No.	Specification(Contents Amount)	Contents Per Box
—	12611	50g	1 Piece

Procare MF Pate

Paste Skin Barrier

Features

- Protects skin from excrement by filling gaps between stoma and faceplate.
- Alcohol-free paste remains soft and reduces skin irritation.



Type	Product Code No.	Specification(Contents Amount)	Contents Per Box
—	12601	100g	1 Piece

Procure Paste

Paste Skin Barrier

Features

- Adhesion of appliance can be improved by applying on skin side of faceplate.
- Laminated tube keeps inside soft and can be used up.




Type	Product Code No.	Specification(Contents Amount)	Contents Per Box
—	14431	50g	1 Piece

Procure Soft Wafer-Ring

Wafer Skin Barrier

Features

- Most suitable to smoothen wrinkles and dents around stoma.
- Soft enough to be molded even with hand.
- Ring-shaped wafer is easy to apply to peristomal skin.
- Double-sided adhesive wafer can be adapted to various ostomy appliances' faceplate.



Type	Product Code No.	Specification(Hole Size)	Contents Per Box	
Ring	20	15591	20mmø	8 Pieces
	30	15592	30mmø	8 Pieces

Youcare Digestive

Procure Digestive

Ilefine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procure Urinary

Cellcare Urinary

Accessory

Procare Soft Wafer-Stick

Wafer Skin Barrier

Features

- Most suitable to smoothen wrinkles and dents around stoma.
- Very soft and flexible enough to be molded freely.
- Double-sided adhesive wafer can be adapted to various ostomy appliances' faceplate.



Type	Product Code No.	Specification(LxW)	Contents Per Box
Stick	15601	15mm×100mm	10 Pieces

Cellcare-Wafer

Wafer Skin Barrier

Features

- Skin barrier wafer contains ceramide, a component that helps preserve moisture in skin.
- Smoothen wrinkles and dents, fitting skin well without collapsing.




Type	Product Code No.	Specification(LxW)	Contents Per Box
10	18391	100mm×100mm	5 Pieces

Procare Wafer

Wafer Skin Barrier

Features

- Most suitable to smoothen wrinkles and dents around stoma.
- Flexibly adjustable to wrinkles and dents by cutting or stacking.



Type	Product Code No.	Specification(LxW)	Contents Per Box
10	12581	100mm×100mm	10 Pieces
20	12582	200mm×195mm	5 Pieces

Trunk Belt

Trunk Belt for Ostomy Appliance

Features

- Prevents ostomy appliance from slipping and falling off by body motion, improving close contact between faceplate and skin.
- Recommended for people who favor activities and who have difficulties in close contact between faceplate and skin.



Type	Product Code No.	Specification		Contents Per Box
		Waist Range	Applicable Products	
1	12691	60~110cm	Procure1-Dc Procure2-Fc Procure1-Uc Procure2-Fc Uro	1 Piece

Youcare Trunk Belt

Trunk Belt for Ostomy Appliance

Features

- Prevents ostomy appliance from slipping and falling off by body motion, improving close contact between faceplate and skin.
- Recommended for people who favor activities and who have difficulties in close contact between faceplate and skin.



Type	Product Code No.	Specification		Contents Per Box
		Waist Range	Applicable Products	
—	15091	60~110cm	Youcare-TDc Cellcare1-TDs Youcare2-TDf Cellcare1-TDc Youcare2-TD Cellcare1-Dc Cap Youcare2-Df Youcare-Uc Youcare2-D Youcare2-U Iifine D-Cap Cellcare1-Uc Iifine TD Youcare2-BC Iifine D(except D70)	1 Piece

Pediatric Trunk Belt

Trunk Belt for Ostomy Appliance

Features

- Can be used with Pediatric Procure 1 only.
- Prevents ostomy appliance from slipping and falling off by body motion, improving close contact between faceplate and skin.
- Easy attach and detach.



Type	Product Code No.	Specification		Contents Per Box
		Waist Range	Applicable Products	
—	14441	30~50cm	Pediatric Procure1-D Pediatric Procure1-U	1 Piece

Youcare Digestive

Procure Digestive

Iifine Digestive

Cellcare Digestive

Others Digestive

Youcare Urinary

Procure Urinary

Cellcare Urinary

Accessory

Youcare
Digestive

Trunk Belt with Faceplate

Trunk Belt for
Ostomy Appliance

Features

- Prevents ostomy appliance from slipping and falling off by body motion, improving close contact between faceplate and skin.
- Recommended for people who favor activities and who have difficulties in close contact between faceplate and skin.



Type	Product Code No.	Specification		Contents Per Box
		Waist Range	Applicable Products	
—	12671	60cm~110cm (Hole Size 63mmø)	Procare1-D Procare1-C Cellcare1-TD(except 70) Cellcare1-D Cap(except 70) Procare1-U Cellcare1-U	1 Set

Procare
Digestive

Ilefine
Digestive

Cellcare
Digestive

Youcare Trunk Belt with Faceplate

Trunk Belt for
Ostomy Appliance

Features

- Prevents ostomy appliance from slipping and falling off by body motion, improving close contact between faceplate and skin.
- Recommended for people who favor activities and who have difficulties in close contact between faceplate and skin.



Type	Product Code No.	Specification		Contents Per Box
		Waist Range	Applicable Products	
—	15081	60cm~110cm (Hole Size 73mmø)	Youcare-TD Youcare-D Youcare-C Ilefine D-Cap Flat(except 60) Pospack-B Youcare-U	1 Set

Others
Digestive

Youcare
Urinary

Procare
Urinary


Cellcare
Urinary

Trunk Belt with Faceplate-II

Trunk Belt for
Ostomy Appliance

Features

- Elastic and soft fixing belt provides comfortableness.
- Improving close contact between faceplate and skin, holding whole faceplate.



*Faceplate is not included.

Type	Product Code No.	Specification		Contents Per Box
		Set Contents	Applicable Products	
—	18061	Youcare Trunk Belt : 1 Piece Fixation Band : 2 Pieces	Procare2-F(except 70) Procare2-FA Cellcare2-F(except LL) Cellcare2-Fc Procare2-FA Uro	1 Set


Product Name	Product Code No.	Specification	Contents Per Box
Fixation Band	18071	—	2 Pieces

Youcare 2-BC

2-piece Ostomy System

Features

- Recommended for bathing or swimming thanks to compact design.
- Durable and easy to wipe off water on surface.
- Flesh-colored cap makes it less noticeable.
- Can be used together with faceplate of Youcare 2.



Type	Product Code No.	Specification(Dome Height)	Contents Per Box
M	16681	20mm	2 Pieces
L	16682	20mm	2 Pieces
65	16683	20mm	2 Pieces

Procare 2-BC

2-piece Ostomy System

Features

- Recommended for bathing or swimming thanks to compact design.
- Durable and easy to wipe off water on surface.
- Flesh-colored cap makes it less noticeable.
- Can be used together with faceplate of Procure 2.



Type	Product Code No.	Specification(Dome Height)	Contents Per Box
30	12181	20mm	5 Pieces
40	12182	20mm	5 Pieces
50	12183	20mm	5 Pieces

Cellcare 2-BC

2-piece Ostomy System

Features

- Recommended for bathing or swimming thanks to compact design.
- Durable and easy to wipe off water on surface.
- Flesh-colored cap makes it less noticeable.
- Can be used together with faceplate of Cellcare 2.



Type	Product Code No.	Specification(Dome Height)	Contents Per Box
M	18291	20mm	2 Pieces
L	18292	20mm	2 Pieces

Youcare
Digestive

Procure
Digestive

Ilefine
Digestive

Cellcare
Digestive

Others
Digestive

Youcare
Urinary

Procure
Urinary

Cellcare
Urinary


Accessory

Mini Close

Bathing Pouch

Features

- Recommended when stool volume is less such as after irrigation or during bathing.
- Adhesive ingredient with less skin irritation is used for faceplate.
- Small and oval faceplate reduces uncomfortableness and does not interfere with body motion.
- Flesh-colored bag makes it less noticeable.




Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size(LxW)	
M	12641	40mmø	~39mm	120mm×110mm	30 Pieces
L	12642	40mmø	~39mm	140mm×130mm	30 Pieces

Mini Pad

Bathing Pad

Features

- Recommended for temporary stoma protection such as during bathing.
- Absorbent pad gently protects stoma.
- Flesh-colored pad makes it less noticeable.



Type	Product Code No.	Specification		Contents Per Box
		Adhesive Part Size	Absorbent Part Size	
—	12651	100mm×100mm	65mm×50mm	30 Pieces

Deofine Powder

Powder Deodorant

Features

- Naturally derived polyphenol prevents odor from generating and diffusing.
- Refreshing mint fragrance makes stoma care more comfortable.
- Easy to carry and convenient thanks to individual packing.



Type	Product Code No.	Specification(Contents Amount)	Contents Per Box
—	17891	1.5g/Piece	30 Pieces

Youcare Legbag

Legbag for Urostomy

Features

- Easily increases urine collection capacity by connecting to stoma bag.
- Reduces burden on application site of ostomy appliance by collecting urine at lower leg.
- Non-woven fabric on skin side absorbs sweat and provides comfortableness.



Type	Product Code No.	Specification (Capacity)	Contents Per Box
—	15371	500mL	5 Pieces

Urinary Bag

Urinary Bag

Features

- Recommended for urine collection at bedside by connecting to stoma bag.
- One-way valve prevents back-flow of urine into stoma bag.
- Easy and one-touch connection mechanism with ALCARE's urinary ostomy appliance.




Type	Product Code No.	Specification			Contents Per Box
		Bag Size (LxW) (Tube Length 1m)	Capacity	Applicable Products	
—	12071	320mmx230mm (Tube Length 1m)	2000mL	Procare1-U Procare1-Uc Procare2-U Pediatric Procare1-U	5 Pieces
1	12072	320mmx230mm (Tube Length 1m)	2000mL	Youcare-U Cellcare1-U Youcare-Uc Cellcare1-Uc Youcare2-U Cellcare2-U	5 Pieces

Ostomy Clip

Ostomy Clip

Features

- Can be used with ALCARE's drainable stoma bag only.



Type	Product Code No.	Specification (Applicable Products)	Contents Per Box
—	90791	Youcare-D Youcare2-Df Youcare2-D Procare1-D Procare1-Dc Procare2-D Ilefine	200 Pieces

Youcare
Digestive

Procare
Digestive

Ilefine
Digestive

Cellcare
Digestive

Others
Digestive

Youcare
Urinary

Procare
Urinary

Cellcare
Urinary

Accessory

Accessory

Youcare
Digestive

Procare
Digestive

Ilefine
Digestive

Cellcare
Digestive

Others
Digestive

Youcare
Urinary

Procare
Urinary

Cellcare
Urinary


Accessory

Connector Tube

Connector Tube

- Can be used with ALCARE's urostomy appliance only.

Features



Type	Product Code No.	Specification (Applicable Products)	Contents Per Box
A	19451	Procare1-U Procare1-Uc Procare2-U Pediatric Procare1-U	10 Pieces
B	19452	Youcare-U Cellcare1-U Youcare-Uc Cellcare1-Uc Youcare2-U Cellcare2-U	10 Pieces

Quick Guide for Accessory

Helpful tips when you want to...	Recommended Products	Features	For			Product Name	Page
			Colostomy	Ileostomy	Urostomy		
change ostomy appliance easily	Remover	Remove ostomy appliance from skin easily	●	●	●	Procare Remover	P25
	Cleanser	Clean peristomal skin easily	●	●	●	Remois Cleanse	P25
protect peristomal skin	Skin Barrier Powder	Protect exposed peristomal skin	●	●	●	Procare Powder	P25
	Skin Barrier Paste		●	●	●	Procare MF Pate	P26
	Skin Barrier Paste	Improve adhesion between faceplate and skin	●	●	●	Procare Paste	P26
	Moldable Skin Barrier	Fill wrinkles and dents around stoma	●	●	●	Procare Soft Wafer-Ring	P26
			●	●	●	Procare Soft Wafer-Stick	P27
	Skin Barrier Wafer		●	●	●	Cellcare-Wafer	P27
●			●	●	Procare Wafer	P27	
make sure ostomy appliance attached firmly	Trunk Belt with Faceplate	Improve close contact between faceplate and skin	●		●	Trunk Belt	P28
			●	●	●	Youcare Trunk Belt	P28
			●		●	Pediatric Trunk Belt	P28
			●	●	●	Trunk Belt with Faceplate	P29
			●	●	●	Youcare Trunk Belt with Faceplate	P29
			●	●	●	Trunk Belt with Faceplate- II	P29
enjoy bath time / hot springs	Bathing Cap	Make ostomy appliance less noticeable	●	●	●	Youcare-BC	P30
			●	●	●	Procare-BC	P30
			●	●	●	Cellcare-BC	P30
	Bathing Pouch		●	●	●	Mini Close	P31
	Bathing Pad		●	●		Mini Pad	P31
prevent odor	Ostomy Deodorant	Reduce stool odor	●	●		Deofine Powder	P31
go out without having concerns about urine	Legbag	Increase urine collection capacity temporarily			●	Youcare Legbag	P32
sleep without concerns	Connecting Bag	Increase urine collection capacity			●	Urinary Bag	P32
have enough accessories you always use	Ostomy Clip	Close outlet of colostomy / ileostomy appliance	●	●		Ostomy Clip	P32
	Connector Tube	Connect urostomy appliance to bag			●	Connector Tube	P33



ALCARE Co., Ltd.

19F, Arcacentral Bldg., 1-2-1 Kinshi, Sumida-ku, Tokyo 130-0013, Japan
TEL.+81-3-5611-7819 FAX.+81-3-5611-7827 MAIL. info@alcare.co.jp
www.alcare.co.jp/en