

# Dermapore Dressing

## Semipermeable Waterproof Wound Dressing

### Features

- Air permeability and waterproofing property protects wound from moisture and contamination.
- Moderate flexibility allows to fit uneven parts.
- Semi-transparent dressing allows observing contamination status of pad without removing.



Type	Product Code No.	Specification		Contents Per Box
		Adhesive Part Size	Absorbent Part Size	
1	11991	4.0 cm × 6.0 cm	2.0 cm × 3.5cm	100 pieces
2	11992	5.0 cm × 7.5 cm	2.0 cm × 4.0cm	50 pieces
3	11993	6.0 cm × 10.0 cm	3.0 cm × 6.5cm	50 pieces
4	11994	10.0 cm × 13.0 cm	4.0 cm × 6.5cm	20 pieces
5	11995	10.0 cm × 16.0 cm	5.0 cm × 10.0cm	20 pieces
6	11996	10.0 cm × 21.0 cm	5.0 cm × 15.0cm	20 pieces
7	11997	10.0 cm × 26.0 cm	5.0 cm × 20.0cm	20 pieces
8	11998	10.0 cm × 31.0 cm	5.0 cm × 25.0cm	20 pieces

EOG Sterilized  
**STERILE**