

Elatex

Elastic Adhesive Bandage

Features

- Tightly twisted 100% cotton allows high elastic and strength compression fixing.
- Rubber adhesive agent achieves strong fixing (Compared to our company).
- Unique adhesive coating technique achieves excellent permeability.
- Latex (natural rubber) is not contained.



Type	Product Code No.	Specification Width × Length (Extension)	Contents Per Pack*
3	11881	2.5 cm × 5 m	12 rolls
4	11882	3.8 cm × 5 m	8 rolls
5	11883	5.0 cm × 5 m	6 rolls
7	11884	7.5 cm × 5 m	4 rolls
10	11885	10.0 cm × 5 m	3 rolls

*Pack : Transparent film shrink pack