


Clearhesive

Hydrocolloid Wound Dressing

Features

- Suture wound dressing using thin hydrocolloid.
- Since it uses a hydrocolloid with excellent water absorption speed, it can be applied immediately after surgery, when exudate and blood are likely to come out.
- It has excellent water resistance and can be used while taking a shower.



Type	Product Code No.	Specification	Contents Per Box
2	19381	5.0 cm × 5.0 cm	30 pieces
3	19382	7.5 cm × 5.0 cm	10 pieces
6	19383	15.0 cm × 5.0 cm	10 pieces
9	19384	22.5 cm × 5.0 cm	10 pieces
12	19385	30.0 cm × 5.0 cm	10 pieces

γ-ray Sterilized
STERILE