

Quickfix

Catheter Fixation Assisted Tape

Features

- By simply placing the catheter on the guide printed on the release paper and wrapping it, you can easily perform firm Ω fixation without putting pressure on the skin.
- The base tape remains as is, and only the catheter tape can be reattached as an integrated type that reduces skin trouble.
- It is pre-cut, so no preparation or processing using scissors is required.



Type	Product Code No.	Specification			Contents Per Box
		Base Tape	Tape For Catheter	Applicable Catheter Size	
1	19391	4.0 cm × 4.0 cm	2.0 cm × 3.8 cm	14 Fr or less	100 pieces
2	19392	4.5 cm × 6.5 cm	2.5 cm × 5.2 cm	20 Fr or less	100 pieces
3	19394	5.0 cm × 10.0 cm	4.0 cm × 7.5 cm	36 Fr or less	100 pieces