

User Review of HN 0.5

Audrey CUGNY
Centre Hospitalier de Libourne

<Challenges>

The Hospital center as only one head support reference for brain and HN treatments.

So, when we want to personalize head patient positions and adapt this patient morphology to the head support used, we have opted for Moldcare solution.

<Solution>

Currently, we are working in an Elekta environment (Linac) and also Elekta contention mask (HeadStep).

We use the **HN-0.5** reference (20x15cm) for our HN treatments only (VMAT technique).

<Conclusion>

The use of Moldcare brings comfort for our patients but it is also an important asset for the quality of the repositioning session after session.

The use of this reference allows us to compensate for the lack of customization possible via different shape of head support.

There are not many choices for head support, so this moldcare makes it possible to compensate for the differences in morphology, so as to optimize the alignment and the curvature of the spine.

Manufacturing during the centering scanner is also a significant profit, very simple and quick for the team.

A.CUGNY

