

Silkytex

Flexible Adhesive Bandage

Silkypore

Flexible Adhesive Gauze

Adhesive tape that combines strength
with low irritancy

- Use of emulsion-type adhesive minimizes chemical irritation
- Two lineups available to cope with various cases

Silkytex
Flexible Adhesive Bandage



Elasticity fitting the flexion site

Silkypore
Flexible Adhesive Gauze



Flexibility and high air permeability gentle to skin

- By introducing “Emulsion Technology”, adhesion of Silkytex and Silkypore will be changed to non-organic solvent acrylic, which will provide further skin friendly adhesive and eco-friendly production process. And it has minimal glue residue and is not prone to peeling due to sweating.

Silkytex

Stretchable for optimal comfort while staying reliably attached.

Fine-tunable from gauze fixation to gentle compression.



Silkytex Flexible Adhesive Bandage

Type	Product Code No.	Specification (W×L)	Contents Per Box
White	1	1.25cm×5m	24 roll
	3	2.5cm×5m	12 roll
	4	3.8cm×5m	8 roll
	5	5.0cm×5m	6 roll
	7	7.5cm×5m	4 roll
	10	10.0cm×5m	3 roll
Beige	3	2.5cm×5m	12 roll
	4	3.8cm×5m	8 roll
	5	5.0cm×5m	6 roll
	7	7.5cm×5m	4 roll
	10	10.0cm×5m	3 roll

Silkypore

Uses highly permeable mesh non-woven fabric that allows skin respiration.

Easily applied to curved areas due to its good stretchability.



Silkypore Flexible Adhesive Gauze

Type	Product Code No.	Specification (W×L)	Contents Per Box
White	1	1.25cm×10m	24 roll
	3	2.5cm×10m	12 roll
	5	5.0cm×10m	6 roll
	7	7.5cm×10m	4 roll
	10	10.0cm×10m	3 roll
	15	15.0cm×10m	2 roll
	30	30.0cm×10m	1 roll
Beige	3	2.5cm×10m	12 roll
	5	5.0cm×10m	6 roll
	7	7.5cm×10m	4 roll
	10	10.0cm×10m	3 roll



ALCARE Co., Ltd.

International Sales Div.

19F, Arcacentral Bldg., 1-2-1 Kinshi, Sumida-ku, Tokyo 130-0013, Japan

TEL.+81-3-5611-7819 FAX.+81-3-5611-7827 MAIL. info@alcare.co.jp

www.alcare.co.jp/e