


Multi Fix

Semipermeable Waterproof Surgical Film



Easy apply,
saves care time

- Easily applied in 3 steps
- Can be applied while wearing surgical gloves
- Adhesive in emulsion form minimizes chemical irritation



The image shows a Multifix Pad.

- Can be attached neatly without wrinkles or gaps with a flexible supporting film.
- Water-proof polyurethane film prevents invasion of water from the outside.



Multi Fix

Semipermeable Waterproof Surgical Film

Type	Product Code No.	Specification	Contents Per Box
1	16386	48mm× 48mm	100 pcs
2	16381	60mm× 70mm	100 pcs
3	16382	80mm×100mm	50 pcs
4	16383	100mm×125mm	50 pcs
5	16384	80mm×180mm	20 pcs
6	16385	100mm×200mm	20 pcs
S-2	16391	150mm×150mm	20 pcs

EOG Sterilized



Multi Fix-Pad

Semipermeable Film Dressing with Absorbent Pad

Type	Product Code No.	Specification		Contents Per Box
		Adhesive Part	Absorptive Part	
2	16921	50mm× 65mm	27mm× 40mm	100 pcs
3	16922	80mm×100mm	40mm× 70mm	20 pcs
4	16923	100mm×125mm	55mm× 75mm	20 pcs
5	16924	80mm×150mm	40mm×110mm	20 pcs
6	16925	100mm×200mm	55mm×150mm	20 pcs
7	16926	100mm×250mm	55mm×200mm	20 pcs

EOG Sterilized



Multi Fix-Roll

Semipermeable Waterproof Film Roll

Type	Product Code No.	Specification	Contents Per Box
5	17821	5.0cm×10m	1 pc
7	17822	7.5cm×10m	1 pc
10	17823	10.0cm×10m	1 pc
15	17824	16.0cm×10m	1 pc



ALCARE Co., Ltd.

International Sales Div.

19F, Arcacentral Bldg., 1-2-1 Kinshi, Sumida-ku, Tokyo 130-0013, Japan

TEL.+81-3-5611-7819 FAX.+81-3-5611-7827 MAIL. info@alcare.co.jp

www.alcare.co.jp/e