

Silkeypore Dressing

Micro-Porous Wound Dressing

Features

- Highly permeable mesh non-woven pad allows uninterrupted breathing of skin.
- Moderate flexibility allows to fit uneven parts.
- Less irritated adhesive allows using for even on delicate skins.



Type	Product Code No.	Specification		Contents Per Box
		Adhesive Part Size	Absorbent Part Size	
1	12008	4.0 cm × 6.0 cm	2.0 cm × 3.5 cm	100 pieces
2	12001	5.0 cm × 7.5 cm	2.0 cm × 4.0 cm	50 pieces
3	12002	6.0 cm × 10.0 cm	3.0 cm × 6.5 cm	50 pieces
4	12003	10.0 cm × 13.0 cm	4.0 cm × 6.5 cm	20 pieces
5	12004	10.0 cm × 16.0 cm	5.0 cm × 10.0 cm	20 pieces
6	12005	10.0 cm × 21.0 cm	5.0 cm × 15.0 cm	20 pieces
7	12006	10.0 cm × 26.0 cm	5.0 cm × 20.0 cm	20 pieces
8	12007	10.0 cm × 31.0 cm	5.0 cm × 25.0 cm	20 pieces

EOG Sterilized
STERILE