

# Multi Fix-Pad

## Semipermeable Film Dressing With Absorbent Pad

### Features

- Flexible support film enables application without being wrinkled, so it prevents from moisture and contamination from outside.
- The absorbent pad covered with non-stick polyethylene net make hard to stick to wound surface.
- Transparent adhesive part enables easy observation of absorbed exudate.



Type	Product Code No.	Specification		Contents Per Box
		Adhesive Part Size	Absorbent Part Size	
2	16921	5.0 cm × 6.5 cm	2.7 cm × 4.0 cm	100 pieces
3	16922	8.0 cm × 10.0 cm	4.0 cm × 7.0 cm	20 pieces
4	16923	10.0 cm × 12.5 cm	5.5 cm × 7.5 cm	20 pieces
5	16924	8.0 cm × 15.0 cm	4.0 cm × 11.0cm	20 pieces
6	16925	10.0 cm × 20.0 cm	5.5 cm × 15.0cm	20 pieces
7	16926	10.0 cm × 25.0 cm	5.5 cm × 20.0cm	20 pieces

EOG Sterilized  
**STERILE**