

# Elatex-S

## Elastic Adhesive Bandage

### Features

- Tightly twisted 100% cotton allows high elastic and strength compression fixing.
- Maintain strong fixing despite of low adhesive strength for reducing peeling irritation.
- Skin-friendly acrylic adhesive allows rarely remaining sticky glue after removal.
- Latex (natural rubber) is not contained.



Type	Product Code No.	Specification Width × Length (Extension)	Contents Per Pack*
3	15331	2.5 cm × 5 m	12 rolls
4	15332	3.8 cm × 5 m	8 rolls
5	15333	5.0 cm × 5 m	6 rolls
7	15334	7.5 cm × 5 m	4 rolls
10	15335	10.0 cm × 5 m	3 rolls

\*Pack : Transparent film shrink pack